

Smart Eye AB

Programmer's Guide

Revision 9.2

Programmer's Guide

Revision 9.2

June 10, 2021

Smart Eye AB

©**Smart Eye AB**

Första Långgatan 28B • 413 27 Gothenburg • Sweden

Phone +46 31 60 61 60

www.smarteye.se

support@smarteye.se

Chapter 1

Abbreviations

SET *Smart Eye Tracker*

SEP *Smart Eye Pro*

WCS World Coordinate System

Contents

| | | |
|----------|--|-----------|
| 1 | Abbreviations | 3 |
| 2 | Introduction | 7 |
| 3 | Smart Eye Data Communication | 8 |
| 3.1 | Smart Eye Packet | 8 |
| 3.1.1 | Packet Specification (TCP/IP, UDP/IP) | 8 |
| 3.1.2 | Data Types | 11 |
| 3.1.3 | SubPacket Ids | 12 |
| 3.2 | SEPacketAPI.h | 12 |
| 3.2.1 | findDataInPacket | 12 |
| 3.2.2 | readValue | 13 |
| 3.2.3 | printPacketContents | 13 |
| 3.2.4 | printSmartEyePacket4 | 13 |
| 3.3 | TCP or UDP? | 14 |
| 3.4 | Real-time or delayed data streams | 14 |
| 3.5 | Logfile specification | 14 |
| 3.6 | A detailed look at some output data | 15 |
| 3.6.1 | Head Position Quality | 15 |
| 3.6.2 | Gaze Direction Quality | 15 |
| 3.6.3 | Gaze Ray, Gaze Origin and Gaze Direction | 16 |
| 3.6.4 | Blinks | 16 |
| 3.6.5 | Saccades and fixations | 18 |
| 3.7 | CAN Output Specification | 18 |
| 3.7.1 | World Intersections | 21 |
| 4 | Smart Eye Remote Control | 22 |
| 4.1 | JSON-RPC | 22 |
| 4.1.1 | JSON-RPC 2.0 specification | 22 |
| 4.1.2 | Connecting to the the JSON-RPC Server | 23 |

| | | |
|----------|---|-----------|
| 4.1.3 | Transportation layer | 23 |
| 4.1.4 | RPC Timings | 24 |
| 4.1.5 | Example DLL | 24 |
| 4.1.6 | JSON-RPC Methods | 24 |
| 4.2 | Exit Codes | 45 |
| 5 | Time Synchronization Specification | 46 |
| 5.1 | Writing a User Time dll | 46 |
| 5.1.1 | startTimeService | 47 |
| 5.1.2 | stopTimeService | 47 |
| 5.1.3 | getCurrentTime | 47 |
| 5.1.4 | timeToString | 48 |
| 6 | Example Applications | 49 |
| 6.1 | Socket Client | 49 |
| 6.2 | Can Client | 49 |
| 6.2.1 | receiveMessage | 49 |
| 6.2.2 | receive<messageType>Message | 49 |
| 6.2.3 | initDriver | 50 |
| 6.3 | TimeService | 50 |
| 6.4 | Python Examples | 50 |
| 6.4.1 | can_client.py | 51 |
| 6.4.2 | socket_client.py | 52 |
| 6.4.3 | client.py | 52 |
| 6.4.4 | external_interfaces.py | 52 |
| 6.4.5 | command_line.py | 53 |
| 6.5 | Broadcast Example | 54 |
| 6.6 | RPC Example | 54 |
| 7 | Scripts | 55 |
| 7.1 | Matlab | 55 |
| 7.1.1 | Log2Matlab | 55 |
| 7.1.2 | Saccade and Fixation Analysis | 55 |
| 7.2 | Scilab | 57 |
| 7.2.1 | Importing Logfiles | 57 |
| A | Smart Eye output data ids | 59 |

| | |
|--|-----------|
| B Smart Eye output data descriptions | 63 |
| C Limitations of Smart Eye Pro Editions | 74 |
| C.1 Smart Eye Tracker | 74 |
| C.1.1 Limitation of JSON-RPCs | 74 |
| C.1.2 Limitation of Output Data | 74 |
| C.1.3 Limitation of World Model | 75 |
| C.2 Smart Sim | 75 |
| C.2.1 Limitation of JSON-RPCs | 75 |
| C.2.2 Limitation of Output Data | 75 |
| C.2.3 Limitation of World Model | 76 |
| C.3 Smart Eye XO | 76 |
| C.3.1 Limitation of JSON-RPCs | 76 |
| C.3.2 Limitation of Output Data | 76 |
| C.3.3 Limitation of World Model | 76 |

Chapter 2

Introduction

This document describes technical details of communication protocols available in *Smart Eye Pro* ([SEP](#)), targeted at integrations between [SEP](#) and other systems. The documentation in this document assumes that the reader is already familiar with programming in C/C++ as example code and datatypes are in most cases exclusively presented in these languages. For general usage and configuration of [SEP](#), please see the *Smart Eye Pro Manual*.



Note

The [SEP](#) protocols will be documented fully but some protocols may be partially unavailable on your system depending on configuration and license used. For details on what limitations apply to specific [SEP](#) editions see [Appendix C, Limitations of Smart Eye Pro Editions](#).

Chapter 3

Smart Eye Data Communication

3.1 Smart Eye Packet

There exist important C/C++ header files that are useful when working with *Smart Eye* network data. The header files are located in the folder `API\include` in the installation path of [SEP](#).

3.1.1 Packet Specification (TCP/IP, UDP/IP)

The basic structure of the data packet produced by the different versions of [SEP](#) is the same. Please see Table [3.1](#) for an illustration of the packet structure.

| | | |
|---------------|------------------|---------------------------------|
| Packet header | | Sync Id SEu32, 4 bytes |
| | | Packet type SEu16, 2 bytes |
| | | Packet length SEu16, 2 bytes |
| Subpacket 1 | Subpacket header | Id SEu16, 2 bytes |
| | | Length SEu16, 2 bytes |
| | Subpacket data | Data |
| Subpacket N | Subpacket header | Id SEu16, 2 bytes |
| | | Length SEu16, 2 bytes |
| | Subpacket data | Data |

Table 3.1: A Schematic View Over the Smart Eye Network Packet Structure

The data packet always starts with the packet header. The packet header declares the type of the packet and the size of the data included in the packet. After the header a number of sub packets follow.

The size (Packet length) declared in the packet header indirectly defines the number of sub packets. The size defines the size in bytes and not the actual number of sub packets. The number of sub packets may differ between two consecutive data packets. The client is responsible for interpreting the data packets correctly.

Each sub packet begins with a sub packet header. The sub packet header contains two pieces of information: id and size (length). The id indicates what data the sub packet contains, for example, SELeftGazeDirection or SEHeadPosition. For a full list of all possible data ids see Table A.1.

Please note that the actual data type (double, int etc.) is not indicated in the packet. The client is expected to look up the data type of the data id. The SocketClient (see Section 6.1) sample that is included with the application illustrates how this may be done.



Note

Please note that data is sent in network byte order. I.e. on an Intel based system the byte order needs to be reversed in order to interpret the data correctly. The sample client also illustrates this.



Note

In very old versions of [SEP](#) the package format was different. Those package types are still supported in the software but not documented here. The package format described here is of type 4.

3.1.2 Data Types

Table 3.2: A List of all the available Smart Eye Data Types.

| DataType | TypeID | Description |
|-----------------------------------|---------------|--|
| SEType_u8 | 0x0000 | Unsigned integer 1 byte |
| SEType_u16 | 0x0001 | Unsigned integer 2 bytes |
| SEType_u32 | 0x0002 | Unsigned integer 4 bytes |
| SEType_s32 | 0x0003 | Signed integer 4 bytes |
| SEType_u64 | 0x0004 | Unsigned integer 8 bytes |
| SEType_f64 | 0x0005 | Floating point 8 bytes |
| SEType_Point2D | 0x0006 | SEType_float x SEType_float y |
| SEType_Vect2D | 0x0007 | SEType_float x SEType_float y |
| SEType_Point3D | 0x0008 | SEType_float x SEType_float y SEType_float z |
| SEType_Vect3D | 0x0009 | SEType_float x SEType_float y SEType_float z |
| SEType_String | 0x000A | SEType_u16 characters (N) SEType_u8 character 1 SEType_u8 character 2 ... SEType_u8 character N |
| SEType_Vector SEType_Vector<T> | 0x000B | SEType_u16 elements (N) SEType_u16 typeld of element 1 Element 1 SEType_u16 typeld of element 2 Element 2 ... SEType_u16 typeld of element N Element N |
| SEType_Struct | 0x000C | SEType_u16 elements (N) SEType_String id of element 1 SEType_u16 typeld of element 1 Element 1 SEType_String id of element 2 SEType_u16 typeld of element 2 Element 2 ... SEType_String id of element N SEType_u16 typeld of element N Element N |
| SEType_WorldIntersection | 0x000D | SEType_u16 intersections (0 1) [SEType_Point3D worldPoint] [SEType_Point3D objectPoint] [SEType_String objectName] |

Table 3.2: A List of all the available Smart Eye Data Types.

| Data Type | TypeID | Description |
|---------------------------|--------|---|
| SEType_WorldIntersections | 0x000E | SEType_u16 intersections (N) [SEType_Point3D worldPoint1] [SEType_Point3D objectPoint1] [SEType_String objectName1] ... [SEType_Point3D worldPointN] [SEType_Point3D objectPointN] [SEType_String objectNameN] |
| SEType_PacketHeader | 0x000F | SEType_u32 syncId SEType_u16 packetType SEType_u16 length |
| SEType_SubPacketHeader | 0x0010 | SEType_u16 id SEType_u16 length |
| SEType_f32 | 0x0011 | Floating point 4 bytes |
| SEType_Matrix3x3 | 0x0012 | A 3x3 Matrix |
| SEType_Matrix2x2 | 0x0013 | A 2x2 Matrix |
| SEType_Quaternion | 0x0014 | SEType_float w SEType_float x SEType_float y SEType_float z |
| SEType_UserMarker | 0x0015 | SEType_u16 exists (0 1) [SEType_s32 error] [SEType_u64 timeStamp] [SEType_u64 cameraClock] [SEType_u8 cameraIdx] [SEType_u64 data] |
| SEType_float | N/A | Alias for either SEType_f32 or SEType_f64 |

3.1.3 SubPacket Ids

Please see the Appendix [A](#) for a list of the possible output data.

3.2 SEPacketAPI.h

3.2.1 findDataInPacket

For all "SEOutputDataIds" there exist a function called "findDataInPacket" that finds the correct subpacket in the body and reads the data from that subpacket. All functions have three parameters except for the functions working with the type "WorldIntersections" (note plural), which have 4 parameters.

```
bool findDataInPacket( const unsigned int& numericalId,
                      const char* pPacket,
                      <SEType>& value )
bool findDataInPacket( const unsigned int& numericalId,
                      const char* pPacket,
                      SEWorldIntersectionStruct value[],
                      SEu16& numberOfIntersections );
```

| Parameter | I/O | Description |
|-------------|-----|---|
| numericalId | I | The SEOutputDataId of the desired data |
| pPacket | I | Pointer to beginning of the packet buffer |
| Value | O | The extracted data, can be any SEType |

| Parameter | I/O | Description |
|-----------------------|-----|--|
| numericalId | I | The SEOutputDataId of SEAllWorldIntersections or SEFilteredAllWorldIntersections |
| pPacket | I | Pointer to beginning of the packet buffer |
| Value | O | The array of extracted SEType_WorldIntersections |
| numberOfIntersections | O | Number of objects that has been extracted |

3.2.2 readValue

For all "SETypes" there exist a function called "readValue" that reads data from a packet header, subpacket header or subpacket body. All functions have three parameters except for the function working with the type WorldIntersections (note plural), which have 4 parameters.

```
void readValue( SEPoint2D& value,
               const char* data,
               int& pos);
void readValue( SEWorldIntersectionStruct value[],
               SEu16& numberOfIntersections,
               const char* data,
               int& pos );
```

| Parameter | I/O | Description |
|-----------|-----|---|
| value | O | The extracted data can be any SEType |
| data | I | Pointer to beginning of the packet buffer |
| pos | I/O | Position in packet buffer |

| Parameter | I/O | Description |
|-----------------------|-----|---|
| value | O | The array of extracted WorldIntersections |
| numberOfIntersections | O | Number of objects that has been extracted |
| data | I | Pointer to beginning of the packet buffer |
| pos | I/O | Position in packet buffer |

3.2.3 printPacketContents

A function that prints all data in a packet to standard out.

```
void printPacketContents( char* pPacket, const SEu16& packetType );
```

| Parameter | I/O | Description |
|------------|-----|---|
| pPacket | I | Pointer to beginning of the packet buffer |
| packetType | I | Type of packet, found in the header |

3.2.4 printSmartEyePacket4

A function that prints all data in a packet of type 4 to standard out.

```
int printSmartEyePacket4( char* pPacket );
```

| Parameter | I/O | Description |
|-----------|-----|---|
| pPacket | I | Pointer to beginning of the packet buffer |

3.3 TCP or UDP?

In a real time context, it is strongly recommended to use UDP over TCP as the UDP socket communication is non blocking. The downside of UDP is that packets are not guaranteed to arrive, but on a local network it is practically guaranteed that they will arrive.

A logger that is not time critical can be implemented in TCP, but still it needs to be ensured that all packets are received and handled as soon as they come in. Failure to handle incoming packets will cause system lag, packet loss or even a system crash. Do not rely on the network buffer! Instead copy the incoming data to your own queues or logs as they come in.



Note

You need to handle all incoming TCP packets as soon as they are received.

Please also keep in mind that the TCP protocol itself may add an unknown and variable amount of lag. The TCP Nagle algorithm ¹ is disabled in the sending end (SEP) causing data packets to be transmitted immediately instead of possibly coalescing them into larger chunks. Still there is no guarantee that this will not happen somewhere else along the way to the receiving end. For this reason TCP is not recommended for real time applications.

3.4 Real-time or delayed data streams

Some of the output values cannot be calculated in real-time and therefore will introduce a delay. If a TCP or UDP stream containing any of these outputs is created, all the outputs will be delayed. The duration of the delay depends on the nature of the specific output value. In the case of blinks for example, the stream will be delayed for about 700 ms since that is about the longest time a natural blink can take. However, the data will still be timestamped at the accurate moment in time.

If all delayed outputs in the data selection list are excluded the data will be sent in real-time. In the UDP/TCP dialog, there is a column that will show if it is real-time or not. In order to get both minimum-lag real-time data as well as outputs from the advanced filters, two separate data streams should be created. One that contains the delayed data and one that contains the real-time data.

Text logs are unaffected by this as they are not real-time logs.

3.5 Logfile specification

The log files generated by SEP are text files with ANSI encoding and tab-separated columns. The first line of the log file is the header containing all the names of the data items and all the following lines contains the values.

For data types with several components (x,y,z) the data item is separated into several columns as

$$HeadPosition.x \rightarrow HeadPosition.y \rightarrow HeadPosition.z.$$

For data types that are a vector it is separated with an appending #i where i is the element position of the vector.

¹ http://en.wikipedia.org/wiki/Nagle's_algorithm

A complex example is the *AllWorldIntersections* data type. It is a vector of *WorldIntersection* which is a struct containing two 3D-points *worldPoint* and *objectPoint* and a string *objectName*. It is written in the log file like this:

```
AllWorldIntersections#0.worldPoint.x → AllWorldIntersections#0.worldPoint.y → AllWorldIntersections#0.worldPoint.z →  
AllWorldIntersections#0.objectPoint.x → AllWorldIntersections#0.objectPoint.y → AllWorldIntersections#0.objectPoint.z →  
AllWorldIntersections#0.objectName → AllWorldIntersections#1.worldPoint.x → AllWorldIntersections#1.worldPoint.y →  
AllWorldIntersections#1.worldPoint.z → AllWorldIntersections#1.objectPoint.x → AllWorldIntersections#1.objectPoint.y →  
AllWorldIntersections#1.objectName → AllWorldIntersections#2.worldPoint.x →  
AllWorldIntersections#2.worldPoint.y → AllWorldIntersections#2.worldPoint.z → AllWorldIntersections#2.objectPoint.x →  
AllWorldIntersections#2.objectPoint.y → AllWorldIntersections#2.objectPoint.z → AllWorldIntersections#2.objectName →
```



Note

Be aware! Some data items can get the value null, this will mean that a double tab will be written to the log file like this "\t\t" (→→). Some tokenizing/splitting functions in different programming languages will skip this empty element and a mismatch in the number of columns compared to the header will occur.



Note

Some values will be written with scientific notations, e.g. 9.52913606973326e-005.



Note

For floating point values a precision of 15 digits is used.

3.6 A detailed look at some output data



Note

It is important to consider the quality value for any output data used. A low quality value may indicate that the associated output data is inaccurate. The quality values are those values whose names end with a Q.

3.6.1 Head Position Quality

The head position quality can take 3 values, a value of 1.0 denotes that the head is found and tracked in at least 2 cameras, while the value 0.5 means that only one camera is tracking. The value 0.0 means that no face is found at all.

Generally speaking, a recommendation is to only use head data with a quality value of 1.0.

3.6.2 Gaze Direction Quality

The GazeDirectionQ value can take any floating point value between 0.0 and 1.0. It depends on the placement of cameras and flashes, the tracking algorithms used and each individuals eye features. This means that it is not easy

to set a general threshold, for example 0.3, where all gaze data below that threshold should be considered unreliable, for all possible setups and subjects.

However, it is too tedious for the user to change this threshold for each specific individual during analysis. So in order to be able to use some generic threshold an attempt is made to normalize this value to the range [0.0, 1.0]. The normalized value 0.0 corresponds to the 1st percentile of all collected quality values of the current tracking session while 1.0 corresponds to the 99th percentile respectively. It takes some samples with good and bad quality measurements before the normalization has built up enough statistics to be effective. Therefore the quality value of the first few seconds of a tracking session may not be comparable to the succeeding quality values.



Note

Normalization is done individually for each camera. This means that the quality levels from the different cameras may differ slightly, especially in the beginning of a tracking session.

The quality value that the system outputs is the one of the camera that is expected to have the best view of the eye. There is no averaging between the cameras. When the system thinks another camera has a better view of the eye, that camera's quality value is the output. Consequently, if the cameras have different normalized quality levels, there might be different quality levels put out of the system.

Normally, this effect is only significant at the very beginning of a tracking session. However, if the system is configured in a way that one of the cameras is used very little, there might be a risk that this camera will provide quality levels in a different range than the other cameras.

A procedure that can be used to avoid this problem is to start each tracking session by looking into each camera once and blink the eyes. This will fill up the normalization buffer with enough statistics.

The gaze direction quality depends on whether the gaze is calculated from pupil or iris detections and uses the quality of the respective detections. In case of pupil-based gaze calculation, this quality is related to how well the detected pupil resembles an ellipse shape and is consequently decreased by disturbances in the pupil edge, such as glints or eyelids.

The iris-based gaze calculation is more complicated. The gaze direction quality reflects the strength of the edge between iris and sclera. This value highly depends on the eye color, illumination and other factors. Therefore the range of this value is more individual.

3.6.3 Gaze Ray, Gaze Origin and Gaze Direction

A "Gaze Ray" is defined by a "Gaze Origin" point in 3D space plus a "Gaze Direction" vector. The "Gaze Direction" vector is a unit vector representing a direction in 3D space.

There is a separate gaze ray for the left and the right eye. These can in turn be combined into a "Cyclops" gaze ray for a virtual eye positioned at the midpoint between the left and the right eye.

Note that no explicit information about the distance to the gaze target is contained in the gaze rays. By estimating the intersection of the left and right gaze rays, some distance information may be recovered, but it will not be very accurate for large distances.

3.6.4 Blinks

The system has a filter that detects blinks by evaluating the measured eyelid opening samples over a period of approximately 700 ms. This means that blinks that last longer than 700 ms will not be considered as blinks.

All samples that belong to a blink will be marked with blink id in the output data. The blink id is an integer that is incremented for each detected blink. Samples that belong to the same blink occurrence will have the same blink id. When no blink is detected, the blink id value is zero. The blink duration may be calculated as the time difference between last and first samples with the same blink id.




Figure 3.1: Gaze Rays for the left (red) and right (green) eye for a top view of a driver looking at a gaze target (symbolised by a yellow star). The 3D (X, Y, Z relative to WCS) Gaze Origin points are drawn as round dots and the unit-length Gaze Direction vectors are drawn as arrows.

As this filter evaluates samples over a period of time, it will cause a delay in the data stream, for which blink data is selected. See section 3.4 on how to get both real time and delayed data.



Note

Blink output creates a lag, so all other data in an output stream that contain blink data will also lag. For a real time data stream this value should be excluded.

The characteristics of a blink can look like the diagram below. Normally the measured eyelid opening does not go all the way down to 0 mm, thus the dashed line just above the time axis.




Figure 3.2: Blink characteristics and measured values

The system measures the characteristics of a blink. The "Closing amplitude" is the difference in eyelid opening before the closing movement and while the eye is closed. The "Opening amplitude" is the difference between the eyelid opening for closed eye and the eyelid opening just after the opening movement. Opening and closing amplitude are not necessarily equal.

"Mid Closing time" is the point in time when the eye is half way down during the closing movement. More precisely it is the time when the eyelid opening has reached 50% of the closing amplitude. To find this time the system fits a typical eyelid movement curve to the measured eyelid opening samples. The exact time is then given by finding the

point on this curve where it reaches 50% of the amplitude. "Mid Opening time" is determined similarly. Both "Mid Closing time" and "Mid Opening time" are reported according to the UserTimeStamp output if a User Time dll has been specified. Otherwise they are reported according to the RealTimeClock output.

"Max closing speed" is the maximum speed of the closing movement, calculated from the fitted curve. Analogously for "Max opening speed".

The values describing blink characteristics are only non-zero at the first sample of each blink occurrence.

While the "blink id" output data is a consensus measurement, that is based on both eyes' eyelid behavior, the values describing the characteristic of the eyelid movement are determined and sent out for each eye separately.

3.6.5 Saccades and fixations

There is also a filter to detect saccades and fixations. The output of these fields are similar to the blink output. The output is an integer that will increment for each saccade or fixation. I.e. it will be a zero when no saccade is occurring and a non-zero integer will indicate that the n:th saccade is occurring. The fixation value will increment after each blink. This output causes a delay since it need to look at future frames to determine if a saccade or a fixation is occurring or not.



Note

Saccades and fixations output creates a lag, so all other data in an output stream will be delayed. For a real time data stream this value needs to be excluded.

For more information on how to adjust the settings for these filters, read the manual for the *Smart Eye* software that you are working with.

3.7 CAN Output Specification

The *Smart Eye* CAN Output interface uses 31 predefined packets containing much of the same information as the data output specification.

All packets have a numerical id (nID in the table below) that is added to a user defined message base for use on the CAN bus. The default message base is 0x400. The message definitions can be found in the `SECANMsgTypes.h` in the *Smart Eye* API as well as in the CANCaseXL-client sample.

The packet data is normally 8 bytes long, except for the `SE_CAN_EYECLOSURE` and `SE_CAN_PUPILDIAMETER` packets which are 4 bytes each.

To save bandwidth all values (except Framenumber and UserTimeStamp) are multiplied by 10 000 before they are sent, which means that metric values will be in 1/10 mm, time measurement in 1/10 ms, and angular values in 1/10 000 radians. Quality values are between 1 and 10000. Most values will be in the form of 16 bit integers, except for Framenumber which has to be a 32 bit integer and the user defined Timestamp which is a 64 bit integer.

On the CAN bus, little endian is used, i.e the least significant byte is always sent first for multiple byte values.

Some messages are sent with a longer delay due to the nature of the data. This group of packages is called "non-realtime messages" and consists of all messages between and including `SE_CAN_SYNCH_NON_REALTIME` to `SE_CAN_RIGHTBLINKOPENINGSPEED`.

**Note**

Output Data is always sent as little endian on the CAN bus.

Table 3.3: CAN Output Data Specification

| Packet ID | nID | Data Byte Type ID | Unit | Description |
|-------------------------------|-----|--|--|--|
| SE_CAN_SYNCH | 1 | [0..3] u32 Framenumber [4..5] u16 Timestamp [6..7] u16 EstimatedDelay | 1/10ms 1/10ms | The time stamp will wrap around at 6.5536 s and the estimated delay will saturate at 6.5535s |
| SE_CAN_USERTIMESTAMP | 2 | [0..7] u64 UserTimeStamp | User Defined | Used to synchronize Smart Eye Pro with external events. |
| SE_CAN_HEADPOSITION | 3 | [0..1] s16 HeadPosition.x [2..3] s16 HeadPosition.y [4..5] s16 HeadPosition.z [6..7] u16 HeadPositionQ | 1/10mm 1/10mm 1/10mm 0-10000 | The head position, scaled down to 1/10 mm. |
| SE_CAN_HEADROTATION | 4 | [0..1] s16 HeadRotation.x [2..3] s16 HeadRotation.y [4..5] s16 HeadRotation.z [6..7] u16 HeadRotationQ | 1/10000 rad 1/10000 rad 1/10000 rad 0-10000 | The head position in rodrigues format, scaled down to 1/10000 radians. |
| SE_CAN_GAZEORIGIN | 5 | [0..1] s16 GazeOrigin.x [2..3] s16 GazeOrigin.y [4..5] s16 GazeOrigin.z [6..7] u16 GazeOriginQ | 1/10mm 1/10mm 1/10mm 0-10000 | The origin of the gaze vector, scaled down to 1/10 mm. |
| SE_CAN_GAZEDIRECTION | 6 | [0..1] s16 GazeDirection.x [2..3] s16 GazeDirection.y [4..5] s16 GazeDirection.z [6..7] u16 GazeDirectionQ | 0-10000 | The unit vector gaze direction, but between -10000 and 10000. |
| SE_CAN_EYECLOSURE | 7 | [0..1] s16 EyeLidOpening [2..3] u16 EyeLidOpeningQ | 1/10mm 0-10000 | The average eye lid opening value. |
| SE_CAN_BOTH EYE CLOSURES | 8 | [0..1] s16 LeftEyeLidOpening [2..3] U16 LeftEyeLidOpeningQ [4..5] s16 RightEyeLidOpening [6..7] u16 RightEyeLidOpeningQ | 1/10mm 0-10000 1/10mm 0-10000 | Both the left and right eye lid opening values with quality values. |
| SE_CAN_PUPILDIAMETER | 9 | [0..1] s16 PupilDiameter [2..3] u16 PupildiameterQ | 1/10mm 0-10000 | The diameter of the pupil in 1/10 mm scale. |
| SE_CAN_BOTH PUPIL DIAMETERS | 10 | [0..1] s16 LeftPupildiameter [2..3] u16 LeftPupildiameterQ [4..5] s16 RightPupildiameter [6..7] u16 RightPupildiameterQ | 1/10mm 0-10000 1/10mm 0-10000 | The diameter of the pupil in 1/10 mm scale for both the left and the right eye. |
| SE_CAN_ASCIIKEY | 11 | [0] u8 ASCIIKeyboardState | | Last keypress value in ASCII format. |
| SE_CAN_WORLD INTERSECTION | 12 | [0..1] u16 ZoneID [2..3] s16 ObjectPoint.x [4..5] s16 ObjectPoint.y [6..7] s16 ObjectPoint.z | 0-65535 1/10mm 1/10mm 1/10mm | Intersection of a World object. The zone ID and the intersected position are reported. |
| SE_CAN_FILTERED GAZEDIRECTION | 13 | [0..1] s16 FilteredGazeDirection.x [2..3] s16 FilteredGazeDirection.y [4..5] s16 FilteredGazeDirection.z | | Filtred gaze direction in rodrigues format and a quality value |

Table 3.3: CAN Output Data Specification

| Packet ID | nID | Data Byte Type ID | Unit | Description |
|---------------------------------------|-----|---|------------------|---|
| | | [6..7] u16 FilteredGazeDirectionQ | 0-10000 | between 0 and 10000. |
| SE_CAN_FILTERED LEFTGAZEDIRECTION | 14 | [0..1] s16 FilteredLeftGazeDirection.x [2..3] s16 FilteredLeftGazeDirection.y [4..5] s16 FilteredLeftGazeDirection.z [6..7] u16 FilteredLeftGazeDirectionQ | 0-10000 | Filtred gaze direction for the left eye in rodrigues format with a quality value. |
| SE_CAN_FILTERED RIGHTGAZEDIRECTION | 15 | [0..1] s16 FilteredRightGazeDirection.x [2..3] s16 FilteredRightGazeDirection.y [4..5] s16 FilteredRightGazeDirection.z [6..7] u16 FilteredRightGazeDirectionQ | 0-10000 | Filtered gaze direction for the right eye in rodrigues format with a quality value. |
| SE_CAN_SYNCH NON_REALTIME | 16 | [0..3] u32 Framenumber [4..5] u16 Timestamp [6..7] u16 EstimatedDelay | 1/10ms 1/10ms | Same as nID 1 but for non-realtime data |
| SE_CAN_USERTIMESTAMP NON_REALTIME | 17 | [0..7] u64 UserTimeStamp | User Defined | Sane as nID 2 but for non-realtime data |
| SE_CAN_BLINK | 18 | [0..3] u32 Blink | | When the blink detection filter detects a blink, this value will be a non-zero integer for the duration of the blink. This integer increases for each blink and can be seen as a unique blink identifier. |
| SE_CAN_SACCADE | 19 | [0..3] u32 Saccade | | When the saccade detection filter detects a saccade, this value will be a non-zero integer for the duration of the saccade. This integer increases for each saccade and can be seen as a unique saccade identifier. |
| SE_CAN_LEFTBLINK CLOSINGMIDTIME | 20 | [0..7] u64 LeftBlinkClosingMidTime | 100ns | See Section 3.6.4. |
| SE_CAN_LEFTBLINK CLOSINGAMPLITUDE | 21 | [0..1] s16 LeftBlinkClosingAmplitude | 1/10mm | See Section 3.6.4. |
| SE_CAN_LEFTBLINK CLOSINGSPEED | 22 | [0..1] s16 LeftBlinkClosingSpeed | 1/10mm/s | See Section 3.6.4. |
| SE_CAN_LEFTBLINK OPENINGMIDTIME | 23 | [0..7] u64 LeftBlinkOpeningMidTime | 100ns | See Section 3.6.4. |
| SE_CAN_LEFTBLINK OPENINGAMPLITUDE | 24 | [0..1] s16 LeftBlinkOpeningAmplitude | 1/10mm | See Section 3.6.4. |
| SE_CAN_LEFTBLINK OPENINGSPEED | 25 | [0..1] s16 LeftBlinkOpeningSpeed | 1/10mm/s | See Section 3.6.4. |
| SE_CAN_RIGHTBLINK CLOSINGMIDTIME | 26 | [0..7] u64 RightBlinkClosingMidTime | 100ns | See Section 3.6.4. |
| SE_CAN_RIGHTBLINK CLOSINGAMPLITUDE | 27 | [0..1] s16 RightBlinkClosingAmplitude | 1/10mm | See Section 3.6.4. |
| SE_CAN_RIGHTBLINK CLOSINGSPEED | 28 | [0..1] s16 RightBlinkClosingSpeed | 1/10mm/s | See Section 3.6.4. |
| SE_CAN_RIGHTBLINK | 29 | [0..7] u64 RightBlinkOpeningMidTime | 100ns | See Section 3.6.4. |

Table 3.3: CAN Output Data Specification

| Packet ID | nID | Data Byte Type ID | Unit | Description |
|---------------------------------------|-----|---------------------------------------|----------|-------------------------------------|
| OPENINGMIDTIME | | | | |
| SE_CAN_RIGHTBLINK OPENINGAMPLITUDE | 30 | [0..1] s16 RightBlinkOpeningAmplitude | 1/10mm | See Section 3.6.4 . |
| SE_CAN_RIGHTBLINK OPENINGSPEED | 31 | [0..1] s16 RightBlinkOpeningSpeed | 1/10mm/s | See Section 3.6.4 . |

For further details refer to the .dbc file in the "Doc/dbc" folder under the installation path.

3.7.1 World Intersections

To get a numerical value as the id of the object the gaze intersects, we use any numerical values in the current closest world intersection object name. To get values for all different intersections, make sure to add a unique numerical value to all Planes, Zones, Spheres and Boxes in the world model. "*Zone222*" will for example give the zone id 222 and "*111Zone333*" will give a zone id of 111333.

The object positions are from the object position value of the Closest World Intersection, but with the unit 1/10 mm.

Chapter 4

Smart Eye Remote Control

The optional module remote control interface allows you to control [SEP](#) from an external client via an JSON-RPC interface. Since this interface is following a open and public specification you can control [SEP](#) from any computer using any OS and programming language of your choice.

4.1 JSON-RPC

4.1.1 JSON-RPC 2.0 specification

Smart Eye follows the JSON-RPC 2.0 specifications which can be found here <http://json-rpc.org/>. A short summary will be given here.

Request

This is the requests that is being sent from the client to the [SEP](#) JSON-RPC server.

```
{
  "jsonrpc" : "2.0",                //REQUIRED
  "method"   : <string>,             //REQUIRED
  "params"   : <json array> or <json object>, //OPTIONAL
  "id"       : <integer> or <string>  //REQUIRED
}
```

Response

To every request a response is given from the server to the client.

```
{
  "jsonrpc" : "2.0",                //REQUIRED
  "result"   : <any json primitive>, //REQUIRED on SUCCESS, MUST NOT EXIST on ERROR
  "error"    : <error object>,       //REQUIRED on ERROR, MUST NOT EXIST on SUCCESS
  "id"       : <integer> or <string>  //REQUIRED
}
```

Notification

A notification is a request that shall NOT get a response.

```
{
  "jsonrpc" : "2.0",           //REQUIRED
  "method"  : <string>,        //REQUIRED
  "params"  : <json array> or <json object>, //OPTIONAL
}
```

Error object

If an error occurs, this is the format of the error object.

```
{
  "code"      : <integer>           //REQUIRED
  "message"   : <string>,           //REQUIRED, SHOULD be one single sentence
  "data"      : <any json primitive> //OPTIONAL, any additional data that server
                                     //may want to attach with the error.
}
```

Pre-defined error codes

These are predefined by the JSON-RPC specification.

| CODE | HEX | MESSAGE | MEANING |
|---------------------|-----------------------------|------------------|---|
| -32700 | 0xFFFF8044 | Parse error | Invalid JSON was received by the server. An error occurred on the server while parsing the JSON text. |
| -32600 | 0xFFFF80A8 | Invalid request | The JSON sent is not a valid Request object. |
| -32601 | 0xFFFF80A7 | Method not found | The method does not exist / is not available. |
| -32602 | 0xFFFF80A6 | Invalid params | Invalid method parameter(s). |
| -32603 | 0xFFFF80A5 | Internal error | Internal JSON-RPC error. |
| -32000 to -32099 | 0xFFFF8300 to 0xFFFF829D | Server error | Reserved for implementation-defined server-errors. |

4.1.2 Connecting to the the JSON-RPC Server

The server is started automatically together with [SEP](#) assuming that the module *Remote Control* is part of the license. It listens for TCP connections on port 8100 by default, but can be changed through the Windows Registry settings.

Each connected client receives a unique id, so it is possible to have multiple connected clients to the Smart Eye JSON-RPC server. The server can also act as a communication hub allowing for connected clients to send notifications between each other.

4.1.3 Transportation layer

The JSON-RPC 2.0 specification says nothing about the transportation layer. That means that the transportation layer is up to the implementor to choose and can therefore vary between different JSON-RPC libraries.

As mentioned in the previous section [SEP](#) is expecting a TCP connection. That is one part of the transportation layer. But since TCP is a streaming protocol there is no sure way to separate between the JSON messages if they are just sent as pure strings. That is why [SEP](#) is expecting the JSON requests and responses to be encoded as a *netstring*.

What netstring encoding does is simply to add the length of the expected string before it together with a ":" and ends the string with a ",", as "5:smart,3:eye,". The full specification can be found here <https://cr.yp.to/proto/netstrings.txt>.

4.1.4 RPC Timings

Sending lots of RPC calls in short order can sometimes cause dead locks in [SEP](#). If you are experiencing this, as a workaround try adding a one second delay between your RPC calls.

4.1.5 Example DLL

One of the ways to communicate with the JSON-RPC server ([SEP](#)) is to use the native-c dll named "SEJsonRpc.dll". The files can be found under `API\example\bin\` within the installation folder of [SEP](#). The communication between client and server is initialized using the function **seInitialize** and is freed using **seFree**. Refer to the header file and to the commands in Section [4.1.6](#) for information about using the dll. A C++ example application using the SEJsonRpc.dll for RPC communication is also available, see Section [6.6](#).

4.1.6 JSON-RPC Methods

The following methods are available on the [SEP](#) JSON-RPC server.

| | | |
|---|---|--|
| ping | setPlaybackSpeedToMax | startCollectSamplesObject |
| getRealTimeClock | setPlaybackSpeedToRealTime | stopCollectSamples |
| getRPCVersion | setPlaybackPosition | clearAllTargetSamples |
| getProductName | setPlaybackStartStopPositions | clearTargetSamples |
| getCameraType | resumePlayback | retrieveTargetStatistics |
| getFirmwareVersions | pausePlayback | retrieveCalibrationResults |
| getState | setPlaybackRepeatOn | calibrateGaze |
| getRecordingState | setPlaybackRepeatOff | applyGazeCalibration |
| startTracking | startCollectPointSamplesAutomatic | clearGazeCalibration |
| stopTracking | stopCollectPointSamplesAutomatic | shutdown |
| startChessboardTracking | clearProfile | keyDown |
| stopChessboardTracking | saveProfile | keyUp |
| setLogSpecification | loadProfile | subscribeToNotification |
| setLogFile | loadWorldModel | unsubscribeToNotification |
| startLog | getWorldModel | sendNotification |
| stopLog | setWorldModel | getCameraImage |
| setRecordingFile | openDataStreamUDP | isCameraConnectedToUSB3 |
| startRecording | closeDataStreamUDP | |
| stopRecording | openDataStreamTCP | |
| setImageSourceCameras | closeDataStreamTCP | |
| setImageSourceRecording | startCollectSamplesWCS | |

Keep alive



ping

Description:

Sends a ping to the server which replies with a pong. This method call will respond as soon as possible, even though the request queue is not empty or there is an ongoing procedure on the server side. This method is useful to call in a keep-alive polling system.

Result:

[pong] <string> The ping response.

Request Example:

```
{"jsonrpc": "2.0", "method": "ping", "id": 0}
```

Response Example:

```
{"jsonrpc": "2.0", "result": "pong", "id": 0}
```

Time synchronization



getRealTimeClock

Description:

Returns the local computer clock represented as a 64-bit value in FILETIME format as an absolute time since January 1, 1601 (UTC).

Result:

[RealTimeClock] <string> Current real time clock in FILETIME format.

Request Example:

```
{"jsonrpc": "2.0", "method": "getRealTimeClock", "id": 0}
```

Response Example:

```
{"jsonrpc": "2.0", "result": "132282210320546339", "id": 0}
```

Server json-rpc interface version

Use this method to assure compability between different versions of client applications and [SEP](#).



getRPCVersion

Description:

Gets the version of the [SEP](#) JSON-RPC interface.

Result:

[major] <int> Major version.
[minor] <int> Minor version.

Request Example:

```
{"jsonrpc": "2.0", "method": "getRPCVersion", "id": 0}
```

Response Example:

```
{"jsonrpc": "2.0", "result": {"major": 1, "minor": 0}, "id": 0}
```



getProductName

Description:

Returns a user friendly product name.

Result:

[productName] <string> The product name.

Request Example:

```
{"jsonrpc": "2.0", "method": "getProductName", "id": 0}
```

Response Example:

```
{"jsonrpc": "2.0", "result": "Smart Eye Pro 7.2", "id": 0}
```



getCameraType

Description:

Returns the live camera type.

Result:

[cameraType] <string> The camera type.

Request Example:

```
{"jsonrpc": "2.0", "method": "getCameraType", "id": 0}
```

Response Example:

```
{"jsonrpc": "2.0", "result": "Aurora", "id": 0}
```



getFirmwareVersions

Description:

Returns an array of firmware versions of the currently connected cameras.

Result:

[firmwareVersions] <string array> Array of firmware versions. **Note:** Availability and response format differs between camera models.

Request Example:

```
{"jsonrpc": "2.0", "method": "getFirmwareVersions", "id": 0}
```

Response Example:

```
{"jsonrpc": "2.0", "result": ["FW: 5.1, FPGA: 5.2", "FW: 5.1, FPGA: 5.2"], "id": 0}
```

Tracker States

The eye tracker can be in different modes/states. The full list of possible states are found in the table below.

| Name | Hex | Description |
|-------------------------|--------|--|
| Idling | 0x0000 | The system is idle, ready to go into any other state. |
| Tracking | 0x0010 | The system is tracking |
| ChessboardTracker | 0x0015 | The system is currently tracking a chessboard |
| CalibratingCameras | 0x0020 | The system is currently performing camera calibration |
| VerifyingCalibration | 0x0021 | The system is verifying the camera calibration |
| CheckingFocus | 0x0030 | The system is adjusting focus and aperture |
| Mtf | 0x0031 | The system is performing a MTF measurements |
| DefiningWCS | 0x0040 | The system is performing manual WCS definition |
| DefiningAutomaticWCS | 0x0041 | The system is in an defining WCS automatically using a chessboard |
| DefiningCameraTiedWCS | 0x0042 | The system is defining wcs manually, but using cameras for position and rotation |
| GazeCalibrating | 0x0050 | The system is performing gaze calibration |
| GazeVerifying | 0x0051 | The system is verifying gaze calibration |
| RecordingProfile | 0x0060 | The system is collecting snapshots for profile creation |
| PlacingMarkers | 0x0061 | The system is placing markers in profile creation |
| ProfileSelectingPoses | 0x0062 | The system is selecting/reviewing/removing poses |
| ProfileReselectingPoses | 0x0063 | The system is reselecting poses |
| RecordingToFile | 0x0070 | [OBSOLETE] The system is recording to disk |
| ScriptedRecordingToFile | 0x0071 | [OBSOLETE] The system is recording to disk (without open UI dialog) |
| DefiningExpectedPose | 0x0080 | The system is defining expected pose |
| WorldMeasurement | 0x0090 | [OBSOLETE] The system is measuring/building a world model |
| AutoDetectLens | 0x00A0 | The system is detecting lenses |



getState

Description:

Gets the [SEP](#) state.

Result:

[state] <int> The current state.

Request Example:

```
{"jsonrpc": "2.0", "method": "getState", "id": 0}
```

Response Example:

```
{"jsonrpc": "2.0", "result": {"state": 4}, "id": 0}
```

Record States

The full list of possible recordings states is found in the table below.

| Name | Hex | Description |
|--------------|--------|--|
| NotRecording | 0x1000 | The system is not currently recording. |
| Recording | 0x1010 | The system is currently recording. |



getRecordingState

Description:

Gets the [SEP](#) recording state.

Result:

[recordingState] <int> The current recording state.

Request Example:

```
{"jsonrpc": "2.0", "method": "getRecordingState", "id": 0}
```

Response Example:

```
{"jsonrpc": "2.0", "result": {"recordingState": 4096}, "id": 0}
```

Tracking



startTracking

Description:

Start tracking mode.

Request Example:

```
{"jsonrpc": "2.0", "method": "startTracking", "id": 0}
```



stopTracking

Description:

Stop tracking mode.

Request Example:

```
{"jsonrpc": "2.0", "method": "stopTracking", "id": 0}
```



startChessboardTracking

Description:

Start the chessboard tracking mode

Request Example:

```
{  
  "jsonrpc" : "2.0",  
  "method"  : "startChessboardTracking",  
  "id"      : 0  
}
```



stopChessboardTracking

Description:

Stop the chessboard tracking mode

Request Example:

```
{
  "jsonrpc" : "2.0",
  "method" : "stopChessboardTracking",
  "id" : 0
}
```

Logging



setLogSpecification

Description:

Sets the output data specification for text logging.

Parameters:

[0] <string> The output data specification, seperated by ';'. See full list of available output data in the Appendix [B.1](#).

Request Example:

```
{"jsonrpc": "2.0", "method": "setLogSpecification", "params": ["ClosestWorldIntersection; GazeDirectionQ"], "id": 0}
```



setLogFile

Description:

Sets the file path of the output log file.

Parameters:

[0] <string> File path to the specified output log file. The path must be reachable by the server computer.

Request Example:

```
{"jsonrpc": "2.0", "method": "setLogFile", "params": ["d:\\logs\\my_log.log"], "id": 0}
```



startLog

Description:

Starts logging a text log.

Note! This command assumes that the log path has already been set by the command **setLogFile**

Request Example:

```
{"jsonrpc": "2.0", "method": "startLog", "id": 0}
```



stopLog

Description:

Stops logging the active text log.

Request Example:

```
{"jsonrpc": "2.0", "method": "stopLog", "id": 0}
```

Recording



setRecordingFile

Description:

Sets the destination to where recorded data will be saved.

Note! If the destination hard drive is not fast enough or full the recording will fail.

Parameters:

[0] <string> File path to the specified output recording file. The path must be reachable by the server computer.

Request Example:

```
{
  "jsonrpc" : "2.0",
  "method"  : "setRecordingFile",
  "params"  : ["d:\\recordings\\my_recording.sma"],
  "id"      : 0
}
```



startRecording

Description:

Starts recording a movie. Returns the remaining recording time (in seconds) calculated as a function of the data rate of the camera system and the available space on the recording disk.

Note! This command assumes that the recording file path has already been set by the command **setRecordingFile**

Parameters:

[0] <int> [OPTIONAL] Compression parameter, can be omitted. Use "0" to turn the compression off and "2" to turn it on (default). The value "1" is a legacy compression type which is no longer supported.

Request Example:

```
{"jsonrpc": "2.0", "method": "startRecording", "id": 0}
```

Request Example:

```
{"jsonrpc": "2.0", "method": "startRecording", "params": [0], "id": 0}
```

Response Example:

```
{"jsonrpc": "2.0", "result": {"remainingRecordingTime": 12}, "id": 0}
```



stopRecording

Description:

Stops recording a movie.

Request Example:

```
{"jsonrpc": "2.0", "method": "stopRecording", "id": 0}
```

Image source



setImageSourceCameras

Description:

Sets live cameras as the image source.

Request Example:

```
{"jsonrpc": "2.0", "method": "setImageSourceCameras", "id": 0}
```



setImageSourceRecording

Description:

Sets a recording as the image source.

Parameters:

[0] <string> File path to the specified recording file that should be played.

Request Example:

```
{
  "jsonrpc" : "2.0",
  "method"  : "setImageSourceRecording",
  "params"  : ["d:\\recordings\\my_recording.sma"],
  "id"      : 0
}
```

Playback



setPlaybackSpeedToMax

Description:

Sets the recording playback speed to max (as fast as the computer can handle).

Request Example:

```
{
  "jsonrpc" : "2.0",
  "method"  : "setPlaybackSpeedToMax",
  "id"      : 0
}
```



setPlaybackSpeedToRealTime

Description:

Sets the recording playback speed to real time (the same frame rate as it was recorded in).

Request Example:

```
{
  "jsonrpc" : "2.0",
  "method" : "setPlaybackSpeedToRealTime",
  "id" : 0
}
```



setPlaybackPosition

Description:

Set position of the playback. Positions starts at 0 for the first frame, and is incremented by one for each subsequent frame in the recording.

Note! The position of a frame is generally not equal to its FrameNumber.

Parameters:

[position] <int> Which zero-indexed position to play from.

Request Example:

```
{
  "jsonrpc" : "2.0",
  "method" : "setPlaybackPosition",
  "params" : [2],
  "id" : 0
}
```



setPlaybackStartStopPositions

Description:

Set start and stop positions of the playback. Positions starts at 0 for the first frame, and is incremented by one for each subsequent frame in the recording.

Note! The position of a frame is generally not equal to its FrameNumber.

Parameters:

[start] <int> Which zero-indexed position the playback should start from

[stop] <int> Which zero-indexed position the playback should stop at

Request Example:

```
{
  "jsonrpc" : "2.0",
  "method" : "setPlaybackStartStopPositions",
  "params" : [1, 2],
  "id" : 0
}
```




resumePlayback

Description:

Resume the playback

Request Example:

```
{
  "jsonrpc" : "2.0",
  "method" : "resumePlayback ",
  "id" : 0
}
```



pausePlayback

Description:

Pause the playback

Request Example:

```
{
  "jsonrpc" : "2.0",
  "method" : "pausePlayback",
  "id" : 0
}
```



setPlaybackRepeatOn

Description:

Set playback to repeat when finished playing

Request Example:

```
{
  "jsonrpc" : "2.0",
  "method" : "setPlaybackRepeatOn",
  "id" : 0
}
```



setPlaybackRepeatOff

Description:

Set playback to not repeat when finished playing

Request Example:

```
{
  "jsonrpc" : "2.0",
  "method" : "setPlaybackRepeatOff",
  "id" : 0
}
```



startCollectPointSamplesAutomatic

Description:

Orders [SEP](#) to start automatic point samples collection

Request Example:

```
{
  "jsonrpc" : "2.0",
  "method" : "startCollectPointSamplesAutomatic",
  "id" : 0
}
```



stopCollectPointSamplesAutomatic

Description:

Orders [SEP](#) to stop automatic point samples collection

Request Example:

```
{
  "jsonrpc" : "2.0",
  "method" : "stopCollectPointSamplesAutomatic",
  "id" : 0
}
```

Profile



clearProfile

Description:

Clears the current profile.

Request Example:

```
{"jsonrpc":"2.0", "method":"clearProfile", "id":0}
```



saveProfile

Description:

Saves a Profile to the specified path.

Parameters:

[0] <string> File path to the specified profile. The path must be reachable by the server computer.

Note! Don't specify the file ending in the path, [SEP](#) will append that automatically. Depending on which profile mode [SEP](#) is running in the file ending will be different.

Request Example:

```
{"jsonrpc":"2.0", "method":"saveProfile", "params":["d:\\profiles\\my_profile"], "id":0}
```



loadProfile

Description:

Load a Profile with the specified path.

Parameters:

[0] <string> File path to the specified profile. The path must be reachable by the server computer.

Request Example:

```
{"jsonrpc": "2.0", "method": "loadProfile", "params": ["d:\\profiles\\my_profile.hm4"], "id": 0}
```

World model



loadWorldModel

Description:

Loads a World Model with the specified path, or the default world model if no path is specified.

Parameters:

[0] <string> [OPTIONAL] File path to a World Model (.sew file). The path must be reachable by the server computer.

Request Example:

```
{"jsonrpc": "2.0", "method": "loadWorldModel", "params": ["d:\\world_models\\my_world_model.sew"], "id": 0}
```



getWorldModel

Description:

Gets the World Model currently loaded in [SEP](#).

Result:

[worldModel] <string> The loaded World Model of [SEP](#) as a string.

Request Example:

```
{"jsonrpc": "2.0", "method": "getWorldModel", "id": 0}
```

Response Example:

```
{
  "jsonrpc" : "2.0",
  "result" : {
    "worldModel" : "Screen : {\n  name = \"MyScreen\"\n  lowerMiddle = 0.0, 0.05, -0.05\n  size = 0.52, 0.326\n  resolution = 1920, 1200\n  calibrationPoints = 4\n}"
  },
  "id" : 0
}
```



setWorldModel

Description:

Sets the World Model in [SEP](#).

Parameters:

[0] <string> A *Smart Eye* World Model.

Request Example:

```
{
  "jsonrpc" : "2.0",
  "method" : "setWorldModel",
  "params" : ["Screen : {\n  name = \"MyScreen\"\n  lowerMiddle = 0.0, 0.05, -0.05\n  size = 0.52, 0.326\n  resolution = 1920, 1200\n  calibrationPoints = 4\n}"],
  "id" : 0
}
```

Communication UDP/TCP



openDataStreamUDP

Description:

Orders [SEP](#) to open up a UDP data stream.

Parameters:

[0] <string> The destination IP adress. If empty "" it will default to the requester's IP adress.
[1] <int> The port.
[2] <string> The output data specification, seperated by ';'. See full list of available output data in the Appendix [B.1](#).

Request Example:

```
{"jsonrpc":"2.0", "method":"openDataStreamUDP", "params":["192.168.0.3", 5100, "GazeDirection;GazeDirectionQ"], "id":0}
```



closeDataStreamUDP

Description:

Orders [SEP](#) to close a specific UDP data stream.

Parameters:

[0] <string> The destination IP adress. If empty "" it will default to the requester's IP adress.
[1] <int> The port.

Request Example:

```
{"jsonrpc":"2.0", "method":"closeDataStreamUDP", "params":["192.168.0.3", 5100], "id":0}
```



openDataStreamTCP

Description:

Orders [SEP](#) to open up a TCP data stream.

Parameters:

- [0] <int> The port.
- [1] <string> The output data specification, seperated by ';'. See full list of available output data in the Appendix [B.1](#).

Request Example:

```
{ "jsonrpc": "2.0", "method": "openDataStreamTCP", "params": [5100, "GazeDirection;GazeDirectionQ"], "id": 0 }
```



closeDataStreamTCP

Description:

Orders [SEP](#) to close a specific TCP data stream.

Parameters:

- [0] <int> The port.

Request Example:

```
{ "jsonrpc": "2.0", "method": "closeDataStreamTCP", "params": [5100], "id": 0 }
```

Gaze Calibration

The methods in this section are used to perform a gaze calibration procedure remotely. A pseudocode gaze calibration procedure could look like this

```
clearAllTargetSamples;
foreach calibration points do
    display calibration point;
    startCollectSamplesWCS;
    sleep for 1000 ms;
    stopCollectSamples;
end
calibrateGaze;
applyGazeCalibration;
foreach calibration points do
    retrieveTargetStatistics;
end
display calibration results;
```



startCollectSamplesWCS

Description:

Orders [SEP](#) to start collecting points for one 3D-point in space given by (x,y,z) position in the World Coordinate System ([WCS](#)). The tracked person is assumed to be looking at that target point for the duration of the collection time.

Parameters:

- [0] <int> Target id. An id to reference the target by when using other calibration requests.
- [1] <double> The x-component of the 3D-point.
- [2] <double> The y-component of the 3D-point.
- [3] <double> The z-component of the 3D-point.
- [4] <int> Collection timeout in ms.

Request Example:

```
{
  "jsonrpc" : "2.0",
  "method" : "startCollectSamplesWCS",
  "params" : [2, 0.143, 0.365, 0.512, 2000],
  "id" : 0
}
```

**startCollectSamplesByTargetName****Description:**

Orders [SEP](#) to start collecting points for a target point with a specific name that is also specified in the loaded World Model. The tracked person is assumed to be looking at that target point for the duration of the collection time.

Parameters:

- [0] <int> Target id. An id to reference the target by when using other calibration requests.
- [1] <string> The target name. The World Model loaded in [SEP](#) must contain a calibration point with this name.
- [4] <int> Collection timeout in ms.

Request Example:

```
{
  "jsonrpc" : "2.0",
  "method" : "startCollectSamplesByTargetName",
  "params" : [2, "calibPoint1", 2000],
  "id" : 0
}
```

**startCollectSamplesObject****Description:**

Orders [SEP](#) to start collecting points for a target point on the surface of a named object in the World Model. The tracked person is assumed to be looking at that target point for the duration of the collection time.

Parameters:

- [0] <int> Target id. An id to reference the target by when using other calibration requests.
- [1] <string> Object name. The object with this name must exist in the World Model loaded in [SEP](#)
- [1] <double> The x-component of the 3D-point in the objects reference system.
- [2] <double> The y-component of the 3D-point in the objects reference system.
- [3] <double> The z-component of the 3D-point in the objects reference system.
- [3] <int> Collection timeout in ms.

Request Example:

```
{
  "jsonrpc" : "2.0",
  "method" : "startCollectSamplesObject",
  "params" : [2, "LeftScreen", 350, 400, 0, 2000],
  "id" : 0
}
```



stopCollectSamples

Description:

Stops the samples collection. Useful when you want to stop collecting points before the specified timeout period.

Request Example:

```
{"jsonrpc": "2.0", "method": "stopCollectSamples", "id": 0}
```



clearAllTargetSamples

Description:

Clears all the samples for all target points.

Request Example:

```
{"jsonrpc": "2.0", "method": "clearAllTargetSamples", "id": 0}
```



clearTargetSamples

Description:

Clears all samples for a single target point.

Parameters:

[0] <int> Target id of the target to be cleared.

Request Example:

```
{"jsonrpc": "2.0", "method": "clearTargetSamples", "params" : [2], "id": 0}
```



retrieveTargetStatistics

Description:

Retrieves statistics for a single target. It is assumed that a gaze calibration has been performed on target calibration point before retrieving the data. Left eye is index 0 in the arrays, right eye is index 1.

Parameters:

[0] <int> Target id.

Result:

| | | |
|------------|----------------|--|
| [targetId] | <int> | Target id |
| [stdDev] | <double array> | The left and right eye standard deviation of the samples around target calibration point. [Degrees] |
| [accuracy] | <double array> | The left and right eye accuracy of the calibration around target calibration point. [Degrees] |
| [errorsxl] | <double array> | The x-component for each sample after calibration relative to the calibration point (left eye). [Degrees] |
| [errorsyl] | <double array> | The y-component for each sample after calibration relative to the calibration point (left eye). [Degrees] |
| [errorsxr] | <double array> | The x-component for each sample after calibration relative to the calibration point (right eye). [Degrees] |
| [errorsyr] | <double array> | The y-component for each sample after calibration relative to the calibration point (right eye). [Degrees] |

Request Example:

```
{"jsonrpc": "2.0", "method": "retrieveTargetStatistics", "params": [2], "id": 0}
```

Response Example:

```
{
  "jsonrpc" : "2.0",
  "result" : {
    "targetId" : 0,
    "stdDev" : [0.45, 0.462],
    "accuracy" : [0.73, 0.7631],
    "errorsxl" : [-0.0621, 0.3357, 0.2444, 0.1469, ... ],
    "errorsyl" : [0.7448, -0.6037, 0.5173, -0.3936, ... ],
    "errorsxr" : [0.7045, 0.3586, 0.3634, 0.4442, ... ],
    "errorsyr" : [0.7086, 0.8151, -0.1517, -0.5735, ... ]
  },
  "id" : 0
}
```



retrieveCalibrationResults

Description:

Retrieve result from gaze calibration. Left eye is index 0 in the arrays, right eye is index 1.

Result:

| | | |
|------------|----------------|--|
| [stdDev] | <double array> | The left and right eyes standard deviation of all collected samples for all target points. [Degrees] |
| [accuracy] | <double array> | The left and right eyes accuracy of all collected samples for all target points. [Degrees] |

Request Example:

```
{
  "jsonrpc" : "2.0",
  "method" : "retrieveCalibrationResults",
  "id" : 0
}
```

Response Example:

```
{"jsonrpc": "2.0", "result": {"stdDev": [0.5769, 0.462], "accuracy": [0.8775, 0.7631]}, "id": 0}
```




calibrateGaze

Description:

Performs gaze calibration calculations based upon the currently collected samples for all target points. This procedure takes several seconds, so it is recommended to wait for the response of this request in an asynchronous manner. Left eye is index 0 in the arrays, right eye is index 1.

Result:

| | | |
|------------|----------------|--|
| [stdDev] | <double array> | The left and right eyes standard deviation of all collected samples for all target points. [Degrees] |
| [accuracy] | <double array> | The left and right eyes accuracy of all collected samples for all target points. [Degrees] |

Request Example:

```
{"jsonrpc": "2.0", "method": "calibrateGaze", "id": 0}
```

Response Example:

```
{
  "jsonrpc" : "2.0",
  "result" : {
    "stdDev" : [0.5769, 0.462],
    "accuracy" : [0.8775, 0.7631]
  },
  "id" : 0
}
```



applyGazeCalibration

Description:

Applies the calculated gaze calibration. The calibration will not be used until this method is called.

Request Example:

```
{"jsonrpc": "2.0", "method": "applyGazeCalibration", "id": 0}
```



clearGazeCalibration

Description:

Clears the current gaze calibration.

Request Example:

```
{"jsonrpc": "2.0", "method": "clearGazeCalibration", "id": 0}
```

Application control



shutdown

Description:

Shuts [SEP](#) down.

Request Example:

```
{"jsonrpc": "2.0", "method": "shutdown", "id": 0}
```



keyDown

Description:

Simulates pushing a key on the keyboard. It will continue to be pushed down until the **keyUp** command is called.

Parameters:

[0] <string> The key value.

Request Example:

```
{"jsonrpc": "2.0", "method": "keyDown", "params": ["a"], "id": 0}
```



keyUp

Description:

Simulates releasing a key on the keyboard.

Parameters:

[0] <string> The key value.

Request Example:

```
{"jsonrpc": "2.0", "method": "keyUp", "params": ["a"], "id": 0}
```

Notifications

[SEP](#) acts as a central notification server. With this functionality you can send generic messages between clients that are connected to [SEP](#).

There is also some built in notifications into [SEP](#) that gets sent everytime something specific happens. Clients can register to these notifications and react upon them. A full list of the internal notifications is found in the table below.

| Notification name | Description |
|-----------------------|---|
| trackingStarted | Tracking started. |
| trackingStopped | Tracking stopped. |
| endOfMovie | Recording playback reached it's end. |
| collectSamplesStarted | Sample collection started. The same notification is used even if it was started with different kinds of collection methods. |
| collectSamplesStopped | Sample collection stopped. |
| calibrationSuccessful | Calibration was successfully calculated. |
| calibrationFailed | Calibration was unsuccessful. |
| calibrationApplied | Calibration was applied. |
| calibrationCleared | Calibration was cleared. |
| recordingStarted | Recording to file started. |
| recordingStopped | Recording to file stopped. |
| recordingError | An error occurred while recording to file. The "params" field of the notification holds additional error information. |



subscribeToNotification

Description:

Subscribes to a certain notification message. [SEP](#) acts as a notification server.

Parameters:

[0] <string> The notification name.

Request Example:

```
{"jsonrpc": "2.0", "method": "subscribeToNotification", "params": ["trackingStarted"], "id": 0}
```



unsubscribeToNotification

Description:

Unsubscribes to a certain notification message.

Parameters:

[0] <string> The notification name.

Request Example:

```
{"jsonrpc": "2.0", "method": "unsubscribeToNotification", "params": ["trackingStarted"], "id": 0}
```



sendNotification

Description:

Sends a notification message to the server ([SEP](#)). The server will then forward the notification message to all clients subscribing to the specific notification name. The forwarded message will be sent as a notification request according to the json-rpc specifications, see [4.1.1](#).

Parameters:

[0] <string> The notification name.

[1] <string> [OPTIONAL] Any notification parameters that should be attached to the message. This will be sent as a single string and will not be processed by the server, the receiver is assumed to know how this string should be interpreted.

Request Example:

```
{
  "jsonrpc" : "2.0",
  "method"  : "sendNotification",
  "params"  : ["stimuliChanged", "computer1;screen2;stim_1.jpg"],
  "id"      : 0
}
```



getCameraImage

Description:

Gets a snapshot from one of the [SEP](#) camera(s). The image is serialized down to a base64 encoded string. **This is not efficient and is not meant to be used as a high frame rate video feed.**

Parameters:

| | | |
|---------------|----------|--|
| [cameraIndex] | <int> | Camera index |
| [scale] | <double> | Scale of the image (must be 0.01 <= scale <= 10.0) |

Result:

| | | |
|---------------|----------|---|
| [cameraIndex] | <int> | Camera index |
| [width] | <int> | Image width |
| [height] | <int> | Image height |
| [stride] | <int> | The stride. (bytes per line) |
| [timeStamp] | <string> | The timestamp (see SETimeStamp) of the image |
| [imageData] | <string> | Raw uint8 image data encoded as base64 string |

Request Example:

```
{
  "jsonrpc" : "2.0",
  "method" : "getCameraImage",
  "params" : [1, 0.01],
  "id" : 0
}
```

Response Example:

```
{
  "jsonrpc" : "2.0",
  "result" : {
    "cameraIndex" : 0,
    "width" : 192,
    "height" : 144,
    "stride" : 192,
    "timeStamp" : "405305513",
    "imageData" : "CgsLCwsLCgoMCwsLCgoLCws ... LDAoKCgsKCgsLCwoLCgsK"
  },
  "id" : 0
}
```

**isCameraConnectedToUSB3****Description:**

Check if a [SEP](#) camera(s) is connected to a USB3 port.

Parameters:

| | | |
|---------------|-------|---------------|
| [cameraIndex] | <int> | Camera index. |
|---------------|-------|---------------|

Result:

| | | |
|-------------------|--------|-------------------------------------|
| [connectedToUSB3] | <bool> | Is camera connected to USB3 or not. |
|-------------------|--------|-------------------------------------|

Request Example:

```
{
  "jsonrpc" : "2.0",
  "method" : "isCameraConnectedToUSB3",
  "params" : [0],
  "id" : 0
}
```

Response Example:

```
{
  "jsonrpc" : "2.0",
  "result" : {
```

```
    "connectedToUSB3" : True  
  },  
  "id" : 0  
}
```

4.2 Exit Codes

The following exit codes may be returned by [SEP](#):

| Code | Details |
|------|---|
| 0 | Normal exit |
| 1000 | Another instance of SEP is already running. |
| 2000 | Smart Eye Tracker start flag is missing. (Smart Eye Tracker (SET) <i>only</i>) |
| 3000 | Camera is not connected. (SET <i>only</i>) |
| 4000 | Camera is not connected to a USB3 port. (SET <i>only</i>) |
| 5000 | Camera is already in use. (SET <i>only</i>) |
| 6000 | Invalid command line parameters. |

Chapter 5

Time Synchronization Specification

The Time Synchronization Module has two main features:

- It enables the user to provide a user defined 64 bits value as a timestamp to each frame by writing a very simple dll. This is called from the *Smart Eye* system at the moment of exposure for each frame. The UserTimestamp is passed on to all logs and makes it easy to synchronize data from the *Smart Eye* system with data from another system utilizing the same timestamp. Instructions on writing the dll follow below.
- The Time Synchronization module also provides information about latency for each frame, which is the time from the image capture to the time the log data leaves the application either on file, TCP, UDP or rs232. This is important in real time applications. Note that this value differs by several 100 microseconds for different data streams, since they are not send at the exact same time.

5.1 Writing a User Time dll

The installation provides an example on how to implement a dll, see API\examples\HighResolutionTime of the [SEP](#) installation directory. The header file specifying the function implemented by the dll is named TimeService.h and found in API\include in the installation directory.

Compiling the HighResolutionTime project results in the default HighResolutionTime.dll, which provides the realtime clock of the PC.

The dll has four functions that may be modified by the user. The `getCurrentTime` function which should return a 64-bit value representing the UserTimestamp. If necessary, start-up and close-down code can be added into the `startTimeService` and `stopTimeService` functions respectively. The `timeToString` should be able to convert the UserTimestamp to a readable value that may be presented in the user interface.

The `getCurrentTime` function of the dll is called upon capture of a video frame from the cameras. The capture is halted while waiting for a response, thus the `getCurrentTime` function must return very fast and must be non-blocking. Ideally only calls to Windows system dlls should be performed and no locks or any other pieces of code that may halt the execution can be utilized.

Also, the `getCurrentTime` function should ideally be reentrant, meaning that more than one thread can run the function simultaneously without hazards. By default the system expects that the function is reentrant, but this may be configured in the Settings→Software...→User Time Stamp→Current Timestamp DLL tab.

The sample code also contains a project called "HighResolutionTimeTest", which is an application for testing the High-ResolutionTime.dll outside [SEP](#). The HighResolutionTime.exe is run om two clients (can be on the same computer or two different), and the usage is:

```
TimeServiceTest.exe <DLLFileName> <offset> <remoteHost> <clientType>
```

- DLLFileName is the full path to the dll.

- offset is an optional time to add to the current time, set this to 0
- remoteHost is the ip address of the other client
- clientType is either a or b, one client is called a and the other b.

At run time [SEP](#) will add an offset to the time elapsed from the midpoint of camera exposure until the time when `getCurrentTime` is called.

5.1.1 `startTimeService`

Start and possibly synchronize the time service, allocate any resources needed. This function is called exactly once (per client) before any other calls are made from that client.

```
TIMESERVICE_API( startTimeService ) ( TTimeServiceHandle * handle );
```

| Parameter | I/O | Description |
|-----------|-----|----------------------------------|
| handle | 0 | Pointer to a time service handle |

5.1.2 `stopTimeService`

Stop time service and release resources. This function is called exactly once (per client) after all other calls have been made.

```
TIMESERVICE_API( stopTimeService ) ( TTimeServiceHandle handle );
```

| Parameter | I/O | Description |
|-----------|-----|---------------------|
| handle | 1 | time service handle |

5.1.3 `getCurrentTime`

Get the current time in user format. A time offset should be added to the result in order to get the actual time of a past or future event. This function can be called from time-critical sections of the code and should avoid potentially blocking calls and return as soon as possible.

```
TIMESERVICE_API( getCurrentTime ) (
    TTimeServiceHandle handle,
    signed __int64 offset,
    TOpaqueUserTime * currentTime);
```

| Parameter | I/O | Description |
|-------------|-----|---|
| handle | 1 | time service handle |
| offset | 1 | offset in 100 ns units |
| currentTime | 0 | current time in user format, adjusted by offset |

5.1.4 timeToString

Convert a user format timestamp to a human readable text format. The text format is intended for presentation and not data logging, so a one-to-one mapping is not required.

```
TIMESERVICE_API(timeToString) (  
    TTimeServiceHandle handle,  
    TOpaqueUserTime currentTime,  
    unsigned long * bufferSize,  
    char * textBuffer);
```

| Parameter | I/O | Description |
|-------------|-----|---|
| handle | I | time service handle |
| currentTime | O | current time in user format, adjusted by offset |
| bufferSize | I/O | size of text buffer in bytes incl null termination (in = max, out = actual) |
| textBuffer | O | time in text format (null terminated) |

Chapter 6

Example Applications

6.1 Socket Client

The SocketClient example implements an application in C++ for connecting to and receiving output data from [SEP](#) via a TCP or UDP socket.

The SocketClient example files are located in the folder `API\Examples\SocketClient` in the installation path of [SEP](#). Please refer to the the README.md file in that folder for further documentation of the SocketClient example project.

6.2 Can Client

All example files are located in the folder `API\Examples\CanClient` in the installation path of [SEP](#).

This client is an example of how to connect using the CANCaseXL box and receiving messages. Parts of this code can be copied into your own client or used this code to build your client from. The program has an infinite loop that

1. Wait for message
2. Check message type
3. Interpret contents of packet(print content and call user function)
4. Return waiting for packets

6.2.1 receiveMessage

This function takes an event from the CANCaseXL box, determines message type and interpret the contents.

```
int receiveMessage( XLevent xlEvent );
```

| Parameter | I/O | Description |
|-----------|-----|---|
| xlEvent | I | A struct containing an event from the CANCaseXL box |

6.2.2 receive<messageType>Message

This function parses the data from a message and prints the contents to standard out. Valid message types are

- Sync
- UserTimeStamp
- HeadPosition
- HeadRotation
- GazeOrigin
- GazeDirection
- EyeClosure
- BothEyeClosure
- PupilDiameter
- BothPupilDiameter

```
int receive<messageType>Message(unsigned char* message);
```

| Parameter | I/O | Description |
|-----------|-----|---|
| message | I | Pointer to the beginning of the message data buffer |

6.2.3 initDriver

A function that initiates the connection to the CANCaseXL box.

```
XLstatus initDriver( CANBusInfo* canInfo, char* appName );
```

| Parameter | I/O | Description |
|-----------|-----|--|
| canInfo | I | A struct containing information about the CAN connection |
| appName | I | Name of the application in the Vector Hardware configuration |

6.3 TimeService

All example files are located in the folder `API\examples\HighResolutionTime` in the installation path of [SEP](#). Read [Chapter 5](#) for more information.

6.4 Python Examples

There is a suite of Python scripts available for usage together with the [SEP](#) application. Below is a list of the available files with a short explanation. The scripts and usage of them are described in details in separate subsections further below. All scripts require that Python 3 is installed. The scripts can be found in chosen installation folder under `API\Examples\PythonExamples`.

| Script | Description |
|------------------------|---|
| can_types.py | Module containing type structure definitions used by can_client.py. |
| can_client.py | Module containing functionality regarding receiving output data on CAN bus. |
| json_log_can | File containing JSON structured CAN output data. |
| socket_types.py | Module containing type structure definitions used by socket_client.py. |
| socket_client.py | Module containing functionality regarding receiving output data on TCP/UDP socket |
| json_log_socket | File containing JSON structured socket output data. |
| client.py | Contains examples on how to use can_client.py and socket_client.py. |
| external_interfaces.py | Module containing functionality regarding communication with SEP over RPC interface |
| command_line.py | Command line interface wrapping functionality in external_interfaces.py. |
| log_data_list.txt | Contains output message names that can be requested over socket connection. |

6.4.1 can_client.py

The can_client module is intended to help setting up a CAN bus connection and start receiving output data. The module can be imported into other scripts and be used according to example below.

Prerequisite:

CAN output enabled in [SEP](#) application according to [SEP](#) manual.

CAN hardware or virtual CAN bus set up to receive CAN output data according to [SEP](#) manual.

Python 3 version installed.

Python package pywin32 installed.

Example on how to use can_client module:

```
import can_client
client = can_client.CanClient()

client.init()
client.receive()
```

6.4.2 socket_client.py

The socket_client module is intended to help setting up a TCP/UDP connection and start receiving output data. The module can be imported into other scripts and be used according to example below.

Prerequisite:

TCP/UDP output enabled in [SEP](#) application according to [SEP](#) manual.

Example on how to use socket_client module:

```
import socket_client
client = socket_client.SocketClient()

client.init("127.0.0.1", 5002, "tcp")
client.connect()
client.receive()
```

6.4.3 client.py

This module contains usage examples of socket_client and can_client modules. It includes the option to log output data between a range of frames along with examples on how to parse output logs json_log_socket and json_log_can.

Prerequisite:

Same as socket_client and can_client modules depending on which mode is specified.

Example on how to use client module:

1. Open the windows command line.
2. Navigate to the folder where the script resides.
3. Run the following to see the options available:

```
py -3 client.py -h
```

4. Specify command line options and run script again. Example:

```
python -3 client.py --mode can
```

6.4.4 external_interfaces.py

The external_interfaces module contains functionality needed to communicate with the [SEP](#) application over RPC interface. It is intended to simplify the process of writing scripts to interact with the [SEP](#) application. Below is a simple example of how to connect and send a ping to the [SEP](#) application. Other commands can be added in the same fashion.

Prerequisite:

[SEP](#) application running.

Example on how to use external_interfaces module:

```
import external_interfaces

IP = "127.0.0.1"
PORT = "8100"

ext_interface = external_interfaces.ExternalInterface()
ext_interface.connect(IP, PORT)
ext_interface.send_ping()
```

6.4.5 command_line.py

The `command_line` module provides a command line interface that can be used to communicate with [SEP](#) application over RPC interface. It wraps the functionality of `external_interfaces.py` and is intended to be used as is. It is also a good reference example on how to use the `external_interfaces` module.

Prerequisite:

[SEP](#) application running.

Example on how to use `command_line` module:

1. Open the windows command line.
2. Navigate to the folder where the script resides.
3. Run from command line with:

```
py -3 command_line.py
```

6.5 Broadcast Example

All example files are located in the folder `API\Examples\BroadcastExample` in the installation path of [SEP](#). Open the project file in Visual C++ and build the executable to test the example. A prebuilt executable is included in the folder `API\Examples\bin` in the installation path of [SEP](#).

This is an example of how to receive broadcasted video feeds from multiple cameras. The program will look up all the Basler GigE cameras connected to the system and go through the necessary steps for receiving a video feed from each camera. The images received will be displayed in a separate window for each camera for a specified number of frames. The port used defaults to the same port value as [SEP](#) but can be specified for each camera.

6.6 RPC Example

The `RPCExample` shows how to use the `SeJsonRPC.dll` (see section [4.1.5 Example DLL](#)) to communicate with Smart Eye Pro from a C++ application. Please refer to the `README.md` file in the `API\Examples\RPCExample` folder for setup and build instructions.

Chapter 7

Scripts

7.1 Matlab

7.1.1 Log2Matlab

Log2Matlab is an application that converts *Smart Eye* log-files into *Matlab*[®]'s binary data-files (.mat). Multiple log files can be entered and the generated mat-files get the same names as the respective log file. Each column in the log file is stored as a matrix with the same name as the column. If the column is a 3D point, 3D vector or an other type with multiple values the matrix has multiple rows. Finally the six rows in an intersection is stored in `worldPoint.x-z` followed by `objectPoint.x-z`.



Note

*If the *Matlab*[®]-file is going to be used in SciLab it is important to set Scilab flag in the settings. This will shorten variable names that are too long since this is not supported by Scilab.*

7.1.2 Saccade and Fixation Analysis

All *Matlab*[®] files are located in the folder `API\utilities\Matlab` in the installation path of [SEP](#).

saccadeAndFixationAnalysis

The file `saccadeAndFixationAnalysis.m` contains a *Matlab*[®] function that analyses the gaze data for saccades and fixations. It draws graphs to show the result of the analysis.

```
SaccadeAndFixationAnalysis( leftGazeDirection, leftGazeDirectionQ, rightGazeDirection,  
rightGazeDirectionQ, timeStamps )
```

| Parameter | I/O | Description |
|---------------------|-----|--|
| leftGazeDirection | I | The leftGazeDirection column from a log. |
| leftGazeDirectionQ | I | The leftGazeDirectionQ column from a log. |
| rightGazeDirection | I | The rightGazeDirection column from a log. |
| rightGazeDirectionQ | I | The rightGazeDirectionQ column from a log. |
| timeStamps | I | The timestamp column from a log. |

toPolarDeg

A *Matlab*[®] function that convert 3D gaze direction vector to polar coordinates in degrees.

```
polarGaze = toPolarDeg(aGazeDirection)
```

| Parameter | I/O | Description |
|----------------|-----|---|
| aGazeDirection | I | A 3D gaze direction vector. |
| polarGaze | O | The 3D gaze direction vector in polar coordinates in degrees. |

removeLowQ

A *Matlab*[®] function that set gaze direction to NaN where the quality of the gaze is less than a threshold.

```
prunedGazeDirection = removeLowQ( gazeDirection, gazeDirectionQ, threshold )
```

| Parameter | I/O | Description |
|---------------------|-----|---|
| gazeDirection | I | A vector of gaze directions. |
| gazeDirectionQ | I | A vector of gaze direction qualities. |
| threshold | I | The threshold for the qualities. |
| prunedGazeDirection | O | A vector of accepted gaze directions and/or NaN values. |

medianFilt

A *Matlab*[®] function that applies a median filter of size length and ignores NaN values that exist in the input data.

```
vecOut = medianFilt(vecIn, length)
```

| Parameter | I/O | Description |
|-----------|-----|------------------------------|
| vecIn | I | A vector of input data. |
| length | I | Size of filter. |
| vecOut | O | A vector of filtered values. |

rotationVelocity

A *Matlab*[®] function that calculates the horizontal (H), vertical (V) and diagonal ($D = \sqrt{H^2 + V^2}$) rotational velocity of a set of gaze directions.

```
vel = rotationVelocity(gazeDirection,timeSec)
```

| Parameter | I/O | Description |
|---------------|-----|------------------------------------|
| gazeDirection | I | A vector of gaze directions. |
| timeSec | I | A vector times in seconds. |
| vel | O | A vector of rotational velocities. |

7.2 Scilab

All *Matlab*[®] files are located in the folder `API\utilities\Scilab` in the installation path of [SEP](#).

7.2.1 Importing Logfiles

The file `getData.sci` contains functions to import data from *Smart Eye* log-files into the program *SciLab*. It is also possible with the module *Log2Matlab* (described in section [7.1.1](#)) to import logs into *Scilab*.

getDataFromFile

Used to load data from a *Smart Eye* log-file into a matrix.

```
[data, headers] = getDataFromFile(filename)
```

| Parameter | I/O | Description |
|-----------|-----|-----------------------------------|
| filename | I | Tab separated smart eye log file. |
| data | O | Matrix of log data. |
| headers | O | Vector of column names. |

getData

Used to load data from a *Smart Eye* log-file, opens a dialog for choosing which log file to open.

```
[data, headers] = getData()
```

| Parameter | I/O | Description |
|-----------|-----|-------------------------|
| data | O | Matrix of log data. |
| headers | O | Vector of column names. |

getColumnByName

Used to get the column vector for the specified name of a column in the log-file, for example `GazeDirectionQ` or `GazeDirection.x`.

```
[col] = getColumnByName(data,header,name)
```

| Parameter | I/O | Description |
|-----------|-----|--------------------------------------|
| data | O | Matrix of <i>Smart Eye Pro</i> data. |
| headers | O | Vector of column names. |
| name | I | The name of the desired log data. |
| col | O | A vector of the desired data. |

getXYZColumnsByName

Used to get the column vectors, one for each component of a 3-dimensional data type, for the specified name, for example `GazeDirection`.

```
[cols] = getXYZColumnsByName(data,header,name)
```

| Parameter | I/O | Description |
|-----------|-----|--|
| data | 0 | Matrix of <i>Smart Eye Pro</i> data. |
| headers | 0 | Vector of column names. |
| name | I | The name of the desired log data. |
| col | 0 | A matrix (3 cols) of the desired data. |

Appendix A

Smart Eye output data ids

Table A.1: A List of all the available Smart Eye Data, their ID and Data Type.

| DataID | EnumValue | DataType |
|--------------------------|-----------|-------------------------------|
| SEFrameNumber | 0x0001 | SEType_u32 |
| SEEstimatedDelay | 0x0002 | SEType_u32 |
| SETimeStamp | 0x0003 | SEType_u64 |
| SEUserTimeStamp | 0x0004 | SEType_u64 |
| SEFrameRate | 0x0005 | SEType_float |
| SECameraPositions | 0x0006 | SEType_Vector<SEType_Point3D> |
| SECameraRotations | 0x0007 | SEType_Vector<SEType_Vect3D> |
| SEUserDefinedData | 0x0008 | SEType_u64 |
| SERealTimeClock | 0x0009 | SEType_u64 |
| SEHeadPosition | 0x0010 | SEType_Point3D |
| SEHeadPositionQ | 0x0011 | SEType_float |
| SEHeadRotationRodrigues | 0x0012 | SEType_Vect3D |
| SEHeadRotationQuaternion | 0x001d | SEType_Quaternion |
| SEHeadLeftEarDirection | 0x0015 | SEType_Vect3D |
| SEHeadUpDirection | 0x0014 | SEType_Vect3D |
| SEHeadNoseDirection | 0x0013 | SEType_Vect3D |
| SEHeadHeading | 0x0016 | SEType_float |
| SEHeadPitch | 0x0017 | SEType_float |
| SEHeadRoll | 0x0018 | SEType_float |
| SEHeadRotationQ | 0x0019 | SEType_float |
| SEGazeOrigin | 0x001a | SEType_Point3D |
| SELeftGazeOrigin | 0x001b | SEType_Point3D |
| SERightGazeOrigin | 0x001c | SEType_Point3D |
| SEEyePosition | 0x0020 | SEType_Point3D |
| SEGazeDirection | 0x0021 | SEType_Vect3D |
| SEGazeDirectionQ | 0x0022 | SEType_float |
| SELeftEyePosition | 0x0023 | SEType_Point3D |
| SELeftGazeDirection | 0x0024 | SEType_Vect3D |
| SELeftGazeDirectionQ | 0x0025 | SEType_float |
| SERightEyePosition | 0x0026 | SEType_Point3D |
| SERightGazeDirection | 0x0027 | SEType_Vect3D |
| SERightGazeDirectionQ | 0x0028 | SEType_float |
| SEGazeHeading | 0x0029 | SEType_float |
| SEGazePitch | 0x002a | SEType_float |
| SELeftGazeHeading | 0x002b | SEType_float |

Table A.1: A List of all the available Smart Eye Data, their ID and Data Type.

| DataID | EnumValue | DataType |
|-------------------------------------|------------------|---------------------------|
| SELeftGazePitch | 0x002c | SEType_float |
| SERightGazeHeading | 0x002d | SEType_float |
| SERightGazePitch | 0x002e | SEType_float |
| SEFilteredGazeDirection | 0x0030 | SEType_Vect3D |
| SEFilteredLeftGazeDirection | 0x0032 | SEType_Vect3D |
| SEFilteredRightGazeDirection | 0x0034 | SEType_Vect3D |
| SEFilteredGazeHeading | 0x0036 | SEType_float |
| SEFilteredGazePitch | 0x0037 | SEType_float |
| SEFilteredLeftGazeHeading | 0x0038 | SEType_float |
| SEFilteredLeftGazePitch | 0x0039 | SEType_float |
| SEFilteredRightGazeHeading | 0x003a | SEType_float |
| SEFilteredRightGazePitch | 0x003b | SEType_float |
| SESaccade | 0x003d | SEType_u32 |
| SEFixation | 0x003e | SEType_u32 |
| SEBlink | 0x003f | SEType_u32 |
| SEClosestWorldIntersection | 0x0040 | SEType_WorldIntersection |
| SEFilteredClosestWorldIntersection | 0x0041 | SEType_WorldIntersection |
| SEAllWorldIntersections | 0x0042 | SEType_WorldIntersections |
| SEFilteredAllWorldIntersections | 0x0043 | SEType_WorldIntersections |
| SEZoneId | 0x0044 | SEType_u16 |
| SEEstimatedClosestWorldIntersection | 0x0045 | SEType_WorldIntersection |
| SEEstimatedAllWorldIntersections | 0x0046 | SEType_WorldIntersections |
| SEHeadClosestWorldIntersection | 0x0049 | SEType_WorldIntersection |
| SEHeadAllWorldIntersections | 0x004a | SEType_WorldIntersections |
| SEEyelidOpening | 0x0050 | SEType_float |
| SEEyelidOpeningQ | 0x0051 | SEType_float |
| SELeftEyelidOpening | 0x0052 | SEType_float |
| SELeftEyelidOpeningQ | 0x0053 | SEType_float |
| SERightEyelidOpening | 0x0054 | SEType_float |
| SERightEyelidOpeningQ | 0x0055 | SEType_float |
| SEKeyboardState | 0x0056 | SEType_String |
| SELeftLowerEyelidExtremePoint | 0x0058 | SEType_Point3D |
| SELeftUpperEyelidExtremePoint | 0x0059 | SEType_Point3D |
| SERightLowerEyelidExtremePoint | 0x005a | SEType_Point3D |
| SERightUpperEyelidExtremePoint | 0x005b | SEType_Point3D |
| SEPupilDiameter | 0x0060 | SEType_float |
| SEPupilDiameterQ | 0x0061 | SEType_float |
| SELeftPupilDiameter | 0x0062 | SEType_float |
| SELeftPupilDiameterQ | 0x0063 | SEType_float |
| SERightPupilDiameter | 0x0064 | SEType_float |
| SERightPupilDiameterQ | 0x0065 | SEType_float |
| SEFilteredPupilDiameter | 0x0066 | SEType_float |
| SEFilteredPupilDiameterQ | 0x0067 | SEType_float |
| SEFilteredLeftPupilDiameter | 0x0068 | SEType_float |
| SEFilteredLeftPupilDiameterQ | 0x0069 | SEType_float |
| SEFilteredRightPupilDiameter | 0x006a | SEType_float |
| SEFilteredRightPupilDiameterQ | 0x006b | SEType_float |
| SEGPSPosition | 0x0070 | SEType_Point2D |
| SEGPSGroundSpeed | 0x0071 | SEType_float |
| SEGPSCourse | 0x0072 | SEType_float |

Table A.1: A List of all the available Smart Eye Data, their ID and Data Type.

| DataID | EnumValue | DataType |
|--|-----------|---------------------------|
| SEGPSTime | 0x0073 | SEType_u64 |
| SEEstimatedGazeOrigin | 0x007a | SEType_Point3D |
| SEEstimatedLeftGazeOrigin | 0x007b | SEType_Point3D |
| SEEstimatedRightGazeOrigin | 0x007c | SEType_Point3D |
| SEEstimatedEyePosition | 0x0080 | SEType_Point3D |
| SEEstimatedGazeDirection | 0x0081 | SEType_Vect3D |
| SEEstimatedGazeDirectionQ | 0x0082 | SEType_float |
| SEEstimatedGazeHeading | 0x0083 | SEType_float |
| SEEstimatedGazePitch | 0x0084 | SEType_float |
| SEEstimatedLeftEyePosition | 0x0085 | SEType_Point3D |
| SEEstimatedLeftGazeDirection | 0x0086 | SEType_Vect3D |
| SEEstimatedLeftGazeDirectionQ | 0x0087 | SEType_float |
| SEEstimatedLeftGazeHeading | 0x0088 | SEType_float |
| SEEstimatedLeftGazePitch | 0x0089 | SEType_float |
| SEEstimatedRightEyePosition | 0x008a | SEType_Point3D |
| SEEstimatedRightGazeDirection | 0x008b | SEType_Vect3D |
| SEEstimatedRightGazeDirectionQ | 0x008c | SEType_float |
| SEEstimatedRightGazeHeading | 0x008d | SEType_float |
| SEEstimatedRightGazePitch | 0x008e | SEType_float |
| SEFilteredEstimatedGazeDirection | 0x0091 | SEType_Vect3D |
| SEFilteredEstimatedGazeDirectionQ | 0x0092 | SEType_float |
| SEFilteredEstimatedGazeHeading | 0x0093 | SEType_float |
| SEFilteredEstimatedGazePitch | 0x0094 | SEType_float |
| SEFilteredEstimatedLeftGazeDirection | 0x0096 | SEType_Vect3D |
| SEFilteredEstimatedLeftGazeDirectionQ | 0x0097 | SEType_float |
| SEFilteredEstimatedLeftGazeHeading | 0x0098 | SEType_float |
| SEFilteredEstimatedLeftGazePitch | 0x0099 | SEType_float |
| SEFilteredEstimatedRightGazeDirection | 0x009b | SEType_Vect3D |
| SEFilteredEstimatedRightGazeDirectionQ | 0x009c | SEType_float |
| SEFilteredEstimatedRightGazeHeading | 0x009d | SEType_float |
| SEFilteredEstimatedRightGazePitch | 0x009e | SEType_float |
| SEASCIKeyboardState | 0x00a4 | SEType_u16 |
| SECalibrationGazeIntersection | 0x00b0 | SEType_WorldIntersection |
| SETaggedGazeIntersection | 0x00b1 | SEType_WorldIntersection |
| SELeftClosestWorldIntersection | 0x00b2 | SEType_WorldIntersection |
| SELeftAllWorldIntersections | 0x00b3 | SEType_WorldIntersections |
| SERightClosestWorldIntersection | 0x00b4 | SEType_WorldIntersection |
| SERightAllWorldIntersections | 0x00b5 | SEType_WorldIntersections |
| SEFilteredLeftClosestWorldIntersection | 0x00b6 | SEType_WorldIntersection |
| SEFilteredLeftAllWorldIntersections | 0x00b7 | SEType_WorldIntersections |
| SEFilteredRightClosestWorldIntersection | 0x00b8 | SEType_WorldIntersection |
| SEFilteredRightAllWorldIntersections | 0x00b9 | SEType_WorldIntersections |
| SEEstimatedLeftClosestWorldIntersection | 0x00ba | SEType_WorldIntersection |
| SEEstimatedLeftAllWorldIntersections | 0x00bb | SEType_WorldIntersections |
| SEEstimatedRightClosestWorldIntersection | 0x00bc | SEType_WorldIntersection |
| SEEstimatedRightAllWorldIntersections | 0x00bd | SEType_WorldIntersections |
| SELeftBlinkClosingMidTime | 0x00e0 | SEType_u64 |
| SELeftBlinkOpeningMidTime | 0x00e1 | SEType_u64 |
| SELeftBlinkClosingAmplitude | 0x00e2 | SEType_float |
| SELeftBlinkOpeningAmplitude | 0x00e3 | SEType_float |

Table A.1: A List of all the available Smart Eye Data, their ID and Data Type.

| DataID | EnumValue | DataType |
|------------------------------|------------------|---------------------------|
| SELeftBlinkClosingSpeed | 0x00e4 | SEType_float |
| SELeftBlinkOpeningSpeed | 0x00e5 | SEType_float |
| SERightBlinkClosingMidTime | 0x00e6 | SEType_u64 |
| SERightBlinkOpeningMidTime | 0x00e7 | SEType_u64 |
| SERightBlinkClosingAmplitude | 0x00e8 | SEType_float |
| SERightBlinkOpeningAmplitude | 0x00e9 | SEType_float |
| SERightBlinkClosingSpeed | 0x00ea | SEType_float |
| SERightBlinkOpeningSpeed | 0x00eb | SEType_float |
| SELeftEyelidState | 0x0390 | SEType_u8 |
| SERightEyelidState | 0x0391 | SEType_u8 |
| SEUserMarker | 0x03a0 | SEType_UserMarker |
| SECameraClocks | 0x03a1 | SEType_Vector<SEType_u64> |

Appendix B

Smart Eye output data descriptions

Table B.1: Descriptions for different Output Data Types

| Name | Unit | Description |
|----------------|-------|---|
| FrameNumber | | <p>A sequential frame number, set by the image capture subsystem. The frame number starts counting from 0 at application start and increments by +1 for each successfully captured image. This means that the frame number of two consecutive frames will always differ by 1 even if the image capture subsystem should happen to drop one or more frames in between.</p> <p>If the frame numbers in the data output should increment by more than +1 from one data record to the next, then it is either an indication of frame loss in the image processing subsystem, probably due to CPU overload or, in the case of UDP communication, loss of network data packets.</p> <p>If an image source error, such as an Ethernet cable temporarily losing connection to its socket, occurs when Smart Eye Pro is running the frame numbering will restart from 0 to indicate that a problem has occurred.</p> <p>N.B: The frame number should not be used to detect frame loss. The only reliable way to determine frame loss is to check the time stamp difference between two consecutive data records.</p> |
| EstimatedDelay | 100ns | <p>The estimated delay from the real-world event (midpoint of image exposure) to the network socket send() system call.</p> <p>[Time Sync]</p> |
| TimeStamp | 100ns | <p>A high-resolution time stamp, measured at the midpoint of image exposure and starting to count from 0 when the application is started. Since this time stamp is based on the PC hardware high-resolution clock, it should only be used to measure time differences over relatively short periods of time.</p> |

Table B.1: Descriptions for different Output Data Types

| Name | Unit | Description |
|------------------------|--------------|---|
| UserTimeStamp | User defined | 64 bits of user defined information, usually an absolute-time time stamp and/or annotation data. The .dll generating this value has to follow the TimeService.h interface (See Chapter 5 for detailed information). [Time Sync] |
| FrameRate | Hz | The current frame rate of the system in frames. Should ideally be 60 or 120 frames per second depending on system configuration. |
| CameraPositions | m | The position of all the cameras defined in the World Coordinate System. |
| CameraRotations | rad | The rodrigues rotation of all the cameras defined in the World Coordinate System. |
| UserDefinedData | | User defined data, not used at the moment. [UserDefinedData] |
| RealTimeClock | 100ns | The real time clock reads the computer clock and represents it as a 64-bit value in FILE-TIME format as an absolute time since January 1, 1601 (UTC). |
| HeadPosition | m | The 3D head position in the defined World Coordinate System. |
| HeadPositionQ | | Quality of the head position [0..1] (0.0 = no head tracking, 0.05..0.1 = Face Detection, 0.2 = Face Refinder, 0.2..1.0 = Head Tracking). |
| HeadRotationRodrigues | rad | The 3D head orientation in Rodrigues format in the defined World Coordinate System. |
| HeadRotationQuaternion | | The 3D head orientation as a quaternion in the defined World Coordinate System. |
| HeadLeftEarDirection | | Unit vector defining the 'left ear' direction. Identical to the x-axis of the head rotation matrix. |
| HeadUpDirection | | Unit vector defining the 'up' direction. Identical to the y-axis of the head rotation matrix. |
| HeadNoseDirection | | Unit vector defining the 'nose' direction. Identical to the z-axis of the head rotation matrix. |
| HeadHeading | rad | The left/right-rotation of the head. Also known as "no"-rotation (See Section on "Definition of Euler Angles" in the SEP manual). |
| HeadPitch | rad | The up/down-rotation of the head. Also known as "yes"-rotation (See Section on "Definition of Euler Angles" in the SEP manual). |
| HeadRoll | rad | The tilt-rotation of the head. Also known as "maybe"-rotation (See Section on "Definition of Euler Angles" in the SEP manual). |
| HeadRotationQ | | Quality of the head orientation [0..1] (The same value as the head position quality). |
| GazeOrigin | m | The consensus of the LeftGazeOrigin and RightGazeOrigin. |
| LeftGazeOrigin | m | The center of the pupil/iris of the left eye, where the gaze vector originates. |
| RightGazeOrigin | m | The center of the pupil/iris of the right eye, where the gaze vector originates. |

Table B.1: Descriptions for different Output Data Types

| Name | Unit | Description |
|---------------------|------|--|
| EyePosition | m | The virtual 3D eye position given in the defined World Coordinate System. |
| GazeDirection | | A unit vector originating in the virtual eye position describing the direction of the gaze. Calculated as the average of the two gaze vectors originating at the physical eyes. If vergence is of interest, the eyes should be handled separately instead. |
| GazeDirectionQ | | Quality of the gaze direction measurement [0..1]. Depends on the quality of the pupil, iris and glint detection and the degree of agreement between measurements from different eye clips. Calculated as the maximum of the quality of the both physical eyes (see Section 3.6.2). |
| LeftEyePosition | m | The 3D position of the center of the left eye ball given in the defined World Coordinate System. |
| LeftGazeDirection | | A unit vector describing the direction of the gaze of the left eye originating from the LeftGazeOrigin. |
| LeftGazeDirectionQ | | Quality of the left gaze direction measurement [0..1]. Depends on the quality of the pupil, iris and glint detection of the left eye and the degree of agreement between measurements from different eye clips for the left eye (see Section 3.6.2). |
| RightEyePosition | m | The 3D position of the center of the right eye ball given in the defined World Coordinate System. |
| RightGazeDirection | | A unit vector describing the direction of the gaze of the right eye, originating from the LeftGazeOrigin. |
| RightGazeDirectionQ | | Quality of the right gaze direction measurement [0..1]. Depends on the quality of the pupil, iris and glint detection of the right eye and the degree of agreement between measurements from different eye clips for the right eye (see Section 3.6.2). |
| GazeHeading | rad | The left/right angle of the GazeDirection (See Section on "Definition of Euler Angles" in the SEP manual). |
| GazePitch | rad | The up/down angle of the GazeDirection (See Section on "Definition of Euler Angles" in the SEP manual). |
| LeftGazeHeading | rad | The left/right angle of the LeftGazeDirection (See Section on "Definition of Euler Angles" in the SEP manual). |
| LeftGazePitch | rad | The up/down angle of the LeftGazeDirection (See Section on "Definition of Euler Angles" in the SEP manual). |
| RightGazeHeading | rad | The left/right angle of the RightGazeDirection (See Section on "Definition of Euler Angles" in the SEP manual). |

Table B.1: Descriptions for different Output Data Types

| Name | Unit | Description |
|-----------------------------------|------|--|
| RightGazePitch | rad | The up/down angle of the RightGazeDirection (See Section on "Definition of Euler Angles" in the SEP manual). |
| Saccade | | When the saccade detection filter detects a saccade, this value will be a non-zero integer for the duration of the saccade. This integer increases for each saccade and can be seen as a unique saccade identifier. |
| Fixation | | When the fixation detection filter detects a fixation, this value will be a non-zero integer for the duration of the fixation. This integer increases for each fixation and can be seen as a unique fixation identifier. |
| Blink | | When the blink detection filter detects a blink, this value will be a non-zero integer for the duration of the blink. This integer increases for each blink and can be seen as a unique blink identifier. |
| ClosestWorldIntersection | | <p>The closest gaze intersection with any of the world objects. The intersection information contains name of object, intersection point in world coordinates and intersection point in object coordinates.</p> <p>Each sub packet of this type starts with an integer indicating the number of world intersections contained in the sub packet. If there are no gaze intersections with world objects for the current frame this integer will be 0, otherwise 1.</p> |
| AllWorldIntersections | | <p>Analogue to ClosestWorldIntersection, but it contains the whole list of intersected world objects. E.g., if the right window of a car is intersected it may be interesting to also find out if the right rear mirror is intersected.</p> <p>As with ClosestWorldIntersections, each sub packet of this type starts with an integer indicating the number of world intersections contained in the sub packet. The difference is that in this case there may be any number of intersections.</p> |
| FilteredAllWorldIntersections | | <p>Analogue to FilteredClosestWorldIntersection, but it contains the whole list of intersected world objects. E.g., if the right window of a car is intersected it may be interesting to also find out if the right rear mirror is intersected. As with ClosestWorldIntersections, each sub packet of this type starts with an integer indicating the number of world intersections contained in the sub packet. The difference is that in this case there may be any number of intersections.</p> |
| NumericalClosestWorldIntersection | | Obsolete. Will be removed. |

Table B.1: Descriptions for different Output Data Types

| Name | Unit | Description |
|-----------------------------------|------|--|
| EstimatedClosestWorldIntersection | | <p>The closest intersection of the estimated gaze vector with any of the world objects. The intersection information contains name of object, intersection point in world coordinates and intersection point in object coordinates.</p> <p>Each sub packet of this type starts with an integer indicating the number of world intersections contained in the sub packet. If there are no gaze intersections with world objects for the current frame this integer will be 0, otherwise 1.</p> |
| EstimatedAllWorldIntersections | | <p>Analogue to EstimatedClosestWorldIntersection, but it contains the whole list of intersected world objects. E.g., if the right window of a car is intersected it may be interesting to also find out if the right rear mirror is intersected.</p> <p>As with ClosestWorldIntersections, each sub packet of this type starts with an integer indicating the number of world intersections contained in the sub packet. The difference is that in this case there may be any number of intersections.</p> |
| HeadClosestWorldIntersection | | <p>The closest intersection of the head nose vector with any of the world objects. The intersection information contains name of object, intersection point in world coordinates and intersection point in object coordinates.</p> <p>Each sub packet of this type starts with an integer indicating the number of world intersections contained in the sub packet. If there are no gaze intersections with world objects for the current frame this integer will be 0, otherwise 1.</p> |
| HeadAllWorldIntersections | | <p>Analogue to HeadClosestWorldIntersection, but it contains the whole list of intersected world objects. E.g., if the right window of a car is intersected it may be interesting to also find out if the right rear mirror is intersected.</p> <p>As with ClosestWorldIntersections, each sub packet of this type starts with an integer indicating the number of world intersections contained in the sub packet. The difference is that in this case there may be any number of intersections.</p> |
| EyelidOpening | m | The average distance between the eyelids of both eyes. |
| EyelidOpeningQ | | <p>Normally in the range 0..1, some subjects may have values larger than 1.0. The value depends on how distinct the eyelids can be detected. The range is individual and can not be compared between subjects. Calculates as the average quality of the both physical eyes.</p> |
| LeftEyelidOpening | m | The distance between the eyelids of the left eye. |

Table B.1: Descriptions for different Output Data Types

| Name | Unit | Description |
|------------------------------|------|--|
| LeftEyelidOpeningQ | | Quality of left eyelid detection (See description for EyelidOpeningQ). |
| RightEyelidOpening | m | The distance between the eyelids of the right eye. |
| RightEyelidOpeningQ | | Quality of right eyelid detection (See description for EyelidOpeningQ). |
| KeyboardState | | The keys that are currently pressed (A-Z, 0-9), useful for marking events in log files or movie recordings. Please note that due to limitations in a normal PC-keyboard, pressing more than two keys at a time gives undefined output. |
| LeftLowerEyelidExtremePoint | m | The midpoint of the eyelid represented in the 3D coordinates of the World Coordinate System. |
| LeftUpperEyelidExtremePoint | m | The midpoint of the eyelid represented in the 3D coordinates of the World Coordinate System. |
| RightLowerEyelidExtremePoint | m | The midpoint of the eyelid represented in the 3D coordinates of the World Coordinate System. |
| RightUpperEyelidExtremePoint | m | The midpoint of the eyelid represented in the 3D coordinates of the World Coordinate System. |
| PupilDiameter | m | The consensus diameter of the pupils of both eyes. [Pupilometry] |
| PupilDiameterQ | | The quality value for the PupilDiameter. The value depends on how distinct the pupil/iris edge can be detected. The value is normalized and will be in the range of [0.0, 1.0]. [Pupilometry] |
| LeftPupilDiameter | m | The diameter of the pupil of the left eye. [Pupilometry] |
| LeftPupilDiameterQ | | The quality value for the LeftPupilDiameter. The value depends on how distinct the pupil/iris edge can be detected. The value is normalized and will be in the range of [0.0, 1.0]. [Pupilometry] |
| RightPupilDiameter | m | The diameter of the pupil of the right eye. [Pupilometry] |
| RightPupilDiameterQ | | The quality value for the RightPupilDiameter. The value depends on how distinct the pupil/iris edge can be detected. The value is normalized and will be in the range of [0.0, 1.0]. [Pupilometry] |
| EstimatedGazeOrigin | m | The calculation of the estimated gaze is based on the estimated eye center as obtained from head tracking. Its origin is the pupil or iris. |

Table B.1: Descriptions for different Output Data Types

| Name | Unit | Description |
|-----------------------------|------|---|
| EstimatedLeftGazeOrigin | m | The calculation of the estimated gaze is based on the estimated eye center as obtained from head tracking. Its origin is the pupil or iris. |
| EstimatedRightGazeOrigin | m | The calculation of the estimated gaze is based on the estimated eye center as obtained from head tracking. Its origin is the pupil or iris. |
| EstimatedEyePosition | m | The position of the consensus eye as obtained from head tracking. |
| EstimatedGazeDirection | | The calculation of the estimated gaze is based on the estimated eye center as obtained from head tracking. This is the direction is the vector from estimated eye center to pupil or iris. |
| EstimatedGazeDirectionQ | | The calculation of the estimated gaze is based on the estimated eye center as obtained from head tracking. The gaze direction quality depends on how distinct the pupil/iris can be detected and how well gaze measurements from different eye clips coincide. |
| EstimatedGazeHeading | rad | The heading angle of the estimated consensus gaze direction. (See Section on "Definition of Euler Angles" in the SEP manual). |
| EstimatedGazePitch | rad | The pitch angle of the estimated consensus gaze direction. (See Section on "Definition of Euler Angles" in the SEP manual). |
| EstimatedLeftEyePosition | m | The position of the left eye as obtained from head tracking. |
| EstimatedLeftGazeDirection | | The calculation of the estimated gaze is based on the estimated eye center as obtained from head tracking. This is the normalized vector from the left estimated eye center to pupil or iris. |
| EstimatedLeftGazeDirectionQ | | The calculation of the estimated gaze is based on the estimated eye center as obtained from head tracking. The gaze direction quality depends on how distinct the pupil/iris can be detected and how well gaze measurements from several left eye clips coincide. |
| EstimatedLeftGazeHeading | rad | The heading angle of the estimated left gaze direction. (See Section on "Definition of Euler Angles" in the SEP manual). |
| EstimatedLeftGazePitch | rad | The pitch angle of the estimated left gaze direction. (See Section on "Definition of Euler Angles" in the SEP manual). |
| EstimatedRightEyePosition | m | The position of the right eye as obtained from head tracking. |
| EstimatedRightGazeDirection | | The calculation of the estimated gaze is based on the estimated eye center as obtained from head tracking. This is the normalized vector from the right estimated eye center to pupil or iris. |

Table B.1: Descriptions for different Output Data Types

| Name | Unit | Description |
|-------------------------------|------|--|
| EstimatedRightGazeDirectionQ | | The calculation of the estimated gaze is based on the estimated eye center as obtained from head tracking. The gaze direction quality depends on how distinct the pupil/iris can be detected and how well gaze measurements from several right eye clips coincide. |
| EstimatedRightGazeHeading | rad | The heading angle of the estimated right gaze direction. (See Section on "Definition of Euler Angles" in the SEP manual). |
| EstimatedRightGazePitch | rad | The pitch angle of the estimated right gaze direction. (See Section on "Definition of Euler Angles" in the SEP manual). |
| ASCIIKeyboardState | | The ASCII code of the pressed key as an integer. |
| CalibrationGazeIntersection | | Only used to send the selected calibration object to clients. Useful if a calibration point shall be displayed on a screen when performing the gaze calibration. |
| TaggedGazeIntersection | | The position of the current gaze target. It contains the name of the object, the position in world coordinates and the position in object coordinates. |
| LeftClosestWorldIntersection | | The closest intersection of the left eye gaze vector with any of the world objects. The intersection information contains name of object, intersection point in world coordinates and intersection point in object coordinates. Each sub packet of this type starts with an integer indicating the number of world intersections contained in the sub packet. If there are no gaze intersections with world objects for the current frame this integer will be 0, otherwise 1. |
| LeftAllWorldIntersections | | Analogue to LeftClosestWorldIntersection, but it contains the whole list of intersected world objects. E.g., if the right window of a car is intersected it may be interesting to also find out if the right rear mirror is intersected. As with LeftClosestWorldIntersection, each sub packet of this type starts with an integer indicating the number of world intersections contained in the sub packet. The difference is that in this case there may be any number of intersections. |
| RightClosestWorldIntersection | | The closest intersection of the right eye gaze vector with any of the world objects. The intersection information contains name of object, intersection point in world coordinates and intersection point in object coordinates. Each sub packet of this type starts with an integer indicating the number of world intersections contained in the sub packet. If there are no gaze intersections with world objects for the current frame this integer will be 0, otherwise 1. |

Table B.1: Descriptions for different Output Data Types

| Name | Unit | Description |
|---------------------------------------|------|--|
| RightAllWorldIntersections | | <p>Analogue to RightClosestWorldIntersection, but it contains the whole list of intersected world objects. E.g., if the right window of a car is intersected it may be interesting to also find out if the right rear mirror is intersected.</p> <p>As with RightClosestWorldIntersection, each sub packet of this type starts with an integer indicating the number of world intersections contained in the sub packet. The difference is that in this case there may be any number of intersections.</p> |
| FilteredLeftClosestWorldIntersection | | <p>The closest intersection of the filtered left gaze vector with any of the world objects. The intersection information contains name of object, intersection point in world coordinates and intersection point in object coordinates.</p> <p>Each sub packet of this type starts with an integer indicating the number of world intersections contained in the sub packet. If there are no gaze intersections with world objects for the current frame this integer will be 0, otherwise 1.</p> |
| FilteredLeftAllWorldIntersections | | <p>Analogue to FilteredLeftClosestWorldIntersection, but it contains the whole list of intersected world objects. E.g., if the right window of a car is intersected it may be interesting to also find out if the right rear mirror is intersected. As with ClosestWorldIntersections, each sub packet of this type starts with an integer indicating the number of world intersections contained in the sub packet. The difference is that in this case there may be any number of intersections.</p> |
| FilteredRightClosestWorldIntersection | | <p>The closest intersection of the filtered right gaze vector with any of the world objects. The intersection information contains name of object, intersection point in world coordinates and intersection point in object coordinates.</p> <p>Each sub packet of this type starts with an integer indicating the number of world intersections contained in the sub packet. If there are no gaze intersections with world objects for the current frame this integer will be 0, otherwise 1.</p> |
| FilteredRightAllWorldIntersections | | <p>Analogue to FilteredRightClosestWorldIntersection, but it contains the whole list of intersected world objects. E.g., if the right window of a car is intersected it may be interesting to also find out if the right rear mirror is intersected. As with ClosestWorldIntersections, each sub packet of this type starts with an integer indicating the number of world intersections contained in the sub packet. The difference is that in this case there may be any number of intersections.</p> |

Table B.1: Descriptions for different Output Data Types

| Name | Unit | Description |
|--|-------|--|
| EstimatedLeftClosestWorldIntersection | | The closest intersection of the estimated left gaze vector with any of the world objects. The intersection information contains name of object, intersection point in world coordinates and intersection point in object coordinates. Each sub packet of this type starts with an integer indicating the number of world intersections contained in the sub packet. If there are no gaze intersections with world objects for the current frame this integer will be 0, otherwise 1. |
| EstimatedLeftAllWorldIntersections | | Analogue to EstimatedLeftClosestWorldIntersection, but it contains the whole list of intersected world objects. E.g., if the right window of a car is intersected it may be interesting to also find out if the right rear mirror is intersected. As with ClosestWorldIntersections, each sub packet of this type starts with an integer indicating the number of world intersections contained in the sub packet. The difference is that in this case there may be any number of intersections. |
| EstimatedRightClosestWorldIntersection | | The closest intersection of the estimated right gaze vector with any of the world objects. The intersection information contains name of object, intersection point in world coordinates and intersection point in object coordinates. Each sub packet of this type starts with an integer indicating the number of world intersections contained in the sub packet. If there are no gaze intersections with world objects for the current frame this integer will be 0, otherwise 1. |
| EstimatedRightAllWorldIntersections | | Analogue to EstimatedRightClosestWorldIntersection, but it contains the whole list of intersected world objects. E.g., if the right window of a car is intersected it may be interesting to also find out if the right rear mirror is intersected. As with ClosestWorldIntersections, each sub packet of this type starts with an integer indicating the number of world intersections contained in the sub packet. The difference is that in this case there may be any number of intersections. |
| LeftBlinkClosingMidTime | 100ns | See Section 3.6.4. |
| LeftBlinkOpeningMidTime | 100ns | See Section 3.6.4. |
| LeftBlinkClosingAmplitude | m | See Section 3.6.4. |
| LeftBlinkOpeningAmplitude | m | See Section 3.6.4. |
| LeftBlinkClosingSpeed | m/s | See Section 3.6.4. |
| LeftBlinkOpeningSpeed | m/s | See Section 3.6.4. |
| RightBlinkClosingMidTime | 100ns | See Section 3.6.4. |
| RightBlinkOpeningMidTime | 100ns | See Section 3.6.4. |
| RightBlinkClosingAmplitude | m | See Section 3.6.4. |

Table B.1: Descriptions for different Output Data Types

| Name | Unit | Description |
|----------------------------|-------|--|
| RightBlinkOpeningAmplitude | m | See Section 3.6.4. |
| RightBlinkClosingSpeed | m/s | See Section 3.6.4. |
| RightBlinkOpeningSpeed | m/s | See Section 3.6.4. |
| LeftEyelidState | | Detected state of left eyelid. Possible values: "Not Set" (0), "Open" (1), or "Closed" (2). |
| RightEyelidState | | Detected state of right eyelid. Possible values: "Not Set" (0), "Open" (1), or "Closed" (2). |
| UserMarker | | User defined marker, containing some data and the specific time it was received. The marker could, for example, represent a key press on some external experiment equipment. The marker also contains an "Error" field, containing an error code from the Smart Eye error specification. The start of each sub-packet of this type contains a 1 if a marker is received, otherwise 0. |
| CameraClocks | 100ns | Camera hardware provided time stamps of image exposure. Characteristics of these time stamps are camera dependent. |

**Note**

The Rodrigues format is a convenient way to describe rotations: The length of the vector defines the rotation angle about the axis defined by the vector it self. Be careful not to interpret the components of the Rodrigues vector as true rotations about the axes in the corresponding coordinate system. This can only be done if all components are close to zero, otherwise not. This vector can easily be converted to a rotation matrix. Example code to do this is shipped with the system.

**Note**

There exist two version of these data items, unfiltered and filtered. The filtered data items have the same name but is preceded with Filtered.

**Note**

There are some regular data items that have estimated values. The estimated data items have the same name but are preceded with Estimated. When using tracking that is not using corneal reflection, the estimated and regular values will display the same value.

Appendix C

Limitations of Smart Eye Pro Editions

C.1 Smart Eye Tracker

C.1.1 Limitation of JSON-RPCs

The following RPC methods are *not* available in [SET](#):

Tracking

startChessboardTracking
stopChessboardTracking

Recording

All methods

Image source

setImageSourceRecording

Playback

All methods

C.1.2 Limitation of Output Data

The following output data ids are *not* available in [SET](#):

SEHeadRotationRodrigues
SEHeadNoseDirection
SEHeadUpDirection
SEHeadLeftEarDirection
SEHeadHeading
SEHeadPitch
SEHeadRoll
SEGazeOrigin
SEGazeHeading
SEGazePitch
SELeftEyePosition
SELeftGazeOrigin
SELeftGazeDirection
SELeftGazeHeading
SELeftGazePitch

SERightEyePosition
SERightGazeOrigin
SERightGazeDirection
SERightGazeHeading
SERightGazePitch
SEFilteredLeftGazeDirection
SEFilteredRightGazeDirection
SEFilteredGazeHeading
SEFilteredGazePitch
SEFilteredLeftGazeHeading
SEFilteredLeftGazePitch
SEFilteredRightGazeHeading
SEFilteredRightGazePitch
SELeftLowerEyelidExtremePoint
SELeftUpperEyelidExtremePoint
SERightLowerEyelidExtremePoint
SERightUpperEyelidExtremePoint
SEEstimatedGazeOrigin
SEEstimatedGazeDirection
SEEstimatedGazeHeading
SEEstimatedGazePitch
SEEstimatedLeftGazeOrigin
SEEstimatedLeftGazeDirection
SEEstimatedLeftGazeHeading
SEEstimatedLeftGazePitch
SEEstimatedRightGazeOrigin
SEEstimatedRightGazeDirection
SEEstimatedRightGazeHeading
SEEstimatedRightGazePitch
SEFilteredEstimatedGazeDirection
SEFilteredEstimatedLeftGazeDirection
SEFilteredEstimatedRightGazeDirection
SEFilteredEstimatedGazeHeading
SEFilteredEstimatedGazePitch
SEFilteredEstimatedLeftGazeHeading
SEFilteredEstimatedLeftGazePitch
SEFilteredEstimatedRightGazeHeading
SEFilteredEstimatedRightGazePitch

C.1.3 Limitation of World Model

No limitation.

C.2 Smart Sim

C.2.1 Limitation of JSON-RPCs

No limitation.

C.2.2 Limitation of Output Data

Smart Sim is limited to *only* the following output data ids:

SEEyePosition

SEHeadPositionQ
SEGazeDirection
SEFilteredGazeDirection
SEGazeDirectionQ
SEFilteredGazeDirectionQ
SEClosestWorldIntersection
SEFilteredClosestWorldIntersection
SEFrameNumber
SERealTimeClock

C.2.3 Limitation of World Model

World model may only contain objects of type Screen.

C.3 Smart Eye XO

C.3.1 Limitation of JSON-RPCs

No limitation.

C.3.2 Limitation of Output Data

No limitation.

C.3.3 Limitation of World Model

No limitation.